



Needle's Eye Snare[®]

TRANSVENOUS RETRIEVAL SET



Capture problematic cardiac leads
for percutaneous extraction.



Needle's Eye Snare[®]

TRANSVENOUS RETRIEVAL SET

The Needle's Eye Snare[®] is a transvenous grasping device that forms a basketlike snare around the targeted cardiac lead or venous catheter body. It is delivered to the vicinity of the lead or catheter through a long, flexible 12 Fr sheath placed coaxially within a larger outer 16 Fr (OD) sheath that has a hemostasis valve at its proximal end.




Encompassing Grasp

Two shape-retaining loops can encircle cardiac leads, indwelling catheters, fragments of catheter tubing or wire guides, and other foreign objects—whether straight or configured and with or without free ends.



Secure Hold



Advancing the inner sheath closes the snare in a controlled manner and ensures a firm grasp of the target object. Inner and outer sheaths may be used to disrupt fibrous binding tissue and, if necessary, support countertraction maneuvers.



Quick Release

Grip can be released by simply retracting the threader mechanism.

Device details



Needle's Eye

This large loop is made of a nitinol alloy that allows it to retain its shape and hook the target object during deployment. Loop profile diameters are available in 13 and 20 mm.

Threader

This smaller nitinol loop passes through the first to fully encircle the target object.

Handle

The design allows the threader to be engaged with one hand.

Plunger

Depressing the plunger advances the threader; retracting it releases the snare.

Side-Arm Flush Port

The side-arm flush port may be used for fluid infusion or to prefill the snare's innermost lumen with saline.

Telescoping 12 Fr PTFE Sheath

Advancing this sheath tightens the snare to firmly grasp the target object.

16 Fr Introducer Sheath

Advancing this sheath over the snared and folded target object prepares it for removal.

Introducer Tube

The introducer tube can slide forward to compress the snare and ease its reinsertion if replacing the telescoping 12 Fr PTFE sheath with the 12 Fr curved femoral sheath.


Luer Connector

Rotating the connector releases the introducer tube from its storage position atop the handle.



Needle's Eye Snare[®]

TRANSVENOUS RETRIEVAL SET



ADDITIONAL FEATURES

- Compatible with the 12 Fr Curved Inner Femoral Sheath for enhanced positioning within the venous system
- Flush port is continuous with the snare's innermost lumen
- Large-diameter Curved Femoral Introducer Sheath (LR-CSS16) comes equipped with Seldinger insertion components
- Highly visible when using fluoroscopy
- Device is supplied sterile and intended for one-time use

Global Product Number	Order Number	Loop Profile Diameter mm	Length cm
Femoral Snare and Workstation			
G26517	LR-NES001	13	94
G26515	LR-NES002	20	94
G26533	LR-NES541	13	54
Femoral Snare			
G26518	LR-SSN001	13	94
G26516	LR-SSN002	20	94
Global Product Number	Order Number	Sheath OD Fr	Sheath Length cm
Components Available Separately for 94 cm Version			
Curved Femoral Introducer Sheath			
G26505	LR-CRVFEM001	12.0	91.4
G26566	LR-CSS16	16.0	68.6

Needle's Eye Snare is a component of Cook Medical's comprehensive approach to percutaneous cardiac lead extraction.



www.cookmedical.com

COOK MEDICAL INCORPORATED
P.O. Box 4195, Bloomington, IN 47402-4195 U.S.A.
Phone: 812 339 2235, Toll Free: 800 457 4500, Toll-Free Fax: 800 554 8335

COOK (CANADA) INC.
111 Sandiford Drive, Stouffville, Ontario, L4A 7X5 CANADA
Phone: 905.640.7110, Toll Free: 800.668.0300

WILLIAM A. COOK AUSTRALIA PTY. LTD.
95 Brandl Street, Brisbane Technology Park, Eight Mile Plains
Brisbane, QLD 4113 AUSTRALIA, Phone: +61 7 3841 1188

WILLIAM COOK EUROPE ApS
Sandet 6, DK-4632, Bjaeverskov, DENMARK, Phone: +45 56 86 86 86

© COOK 2010 LM-DM-NESBR-EN-201009

AORTIC INTERVENTION

CRITICAL CARE

ENDOSCOPY

INTERVENTIONAL RADIOLOGY

LEAD MANAGEMENT

PERIPHERAL INTERVENTION

SURGERY

UROLOGY

WOMEN'S HEALTH