



## PRODUCTS FOR BEDSIDE PROCEDURES





# Contents

## Vascular Access

- 6-7 Spectrum® Minocycline/Rifampin Catheter

## Tracheostomy

- 8-9 Ciaglia Blue Rhino® G2 Advanced Percutaneous Tracheostomy Introducer
- 10-11 Dolphin BT™ Ciaglia Balloon-Assisted Tracheostomy Introducer

## Self-Advancing Enteral Feeding

- 12-13 Tiger 2™ Self-Advancing Nasal Jejunal Feeding Tube

## Difficult Airway

- 14-15 Airway Exchange Catheters
- 16 Frova Intubating Introducer  
Aintree Intubation Catheter
- 17 Melker Emergency Cricothyrotomy Catheter  
Enk Oxygen Flow Modulator

## Centesis & Drainage

- 18-19 Thal-Quick Chest Tube Introducer
- 20 Wayne Pneumothorax Catheter  
Lock Pericardiocentesis Catheter
- 21 Fuhrman Pleural/Pneumopericardial Drainage Catheter  
Peritoneal Lavage Catheter

## **Recognize the benefits of a bedside procedure.**

### **Safety**

Adverse effects such as arrhythmia, cardiac arrest and hypoxaemia may occur in up to 70% of transports.<sup>1</sup> Performing procedures at the bedside improves quality of care through the prevention of safety risks to critically ill patients during movement to the operating room and other departments within the hospital.

### **Efficiency**

Bedside procedures also eliminate lag time spent waiting for an OR and enable more efficient utilization of hospital resources and staff. This may result in quicker treatment administration and an improved ICU bed turnover rate.

### **Cost Savings**

The average U.S. operating room fee is \$66/minute, with rates continuously on the rise.<sup>2</sup> Bedside procedures minimize overall expenses for the hospital by avoiding high OR-related costs while maintaining clinician reimbursement levels.






# Spectrum<sup>®</sup>

MINOCYCLINE/RIFAMPIN CATHETER

## Proven lowest infection rates.

More than two decades of evidence, including over 21 peer-reviewed studies and meta-analyses, prove Spectrum's ability to help prevent catheter-related bloodstream infections (CRBSIs), without causing bacterial resistance. That's evidence of protection no other process or technology can match. As an additional benefit, Spectrum CVCs are available in comprehensive procedural trays with all the necessary components for successful and efficient bedside insertion.

First Trial of Second Generation AGB+<sup>®</sup> and Spectrum



## Proven cost savings.

The cost to treat just one line infection can range from \$11,971<sup>5</sup> to \$56,000.<sup>6</sup> And with 250,000 cases of CRBSI reported in the U.S. each year,<sup>6</sup> these infections can create a significant economic burden for hospitals. In a challenging clinical environment, a hospital that switches to Spectrum may expect to see a decrease in CRBSI rates, attributable mortality and CRBSI-related costs. What might not be so apparent is that even high-performing hospitals can switch to Spectrum to drive incremental improvement in CRBSI rates and still achieve substantial reductions in mortality and cost.



Order Number	Reference Part Number	Catheter Fr	Equivalent gage	Catheter Length cm	Wire Guide Diameter inch	Fits Needle gage
<b>7.0 Fr Triple Lumen</b>						
Curved Needle with Suture and Holder						
G49803	C-UTLMY-701J-ABRM-HC-IHI-FST-A-RD	7.0	13	15	.032	18 TW
G49804	C-UTLMY-701J-RSC-ABRM-HC-IHI-FST-A-RD	7.0	13	20	.032	18 TW
G49807	C-UTLMY-701J-LSC-ABRM-HC-IHI-FST-A-RD	7.0	13	25	.032	18 TW
<b>8.0 Fr Double Lumen</b>						
Curved Needle with Suture and Holder						
G50821	C-UDLMY-801J-ABRM-HC-IHI-FST-A	8.0	12	15	.032	18 TW
G50822	C-UDLMY-801J-RSC-ABRM-HC-IHI-FST-A	8.0	12	20	.032	18 TW
<b>9.0 Fr Triple Lumen</b>						
Curved Needle with Suture and Holder						
G51589	C-UTLMY-901J-ABRM-HC-IHI-FST-A-RD	9.0	11	15	.032	18 TW
G51590	C-UTLMY-901J-RSC-ABRM-HC-IHI-FST-A-RD	9.0	11	20	.032	18 TW
G51591	C-UTLMY-901J-LSC-ABRM-HC-IHI-FST-A-RD	9.0	11	25	.032	18 TW
<b>10.0 Fr Five Lumen</b>						
Curved Needle with Suture and Holder						
G51546	C-UQLMY-1001J-ABRM-HC-IHI-FST-A-RD	10.0	10	15	.032	18 TW
G51548	C-UQLMY-1001J-RSC-ABRM-HC-IHI-FST-A-RD	10.0	10	20	.032	18 TW
G51547	C-UQLMY-1001J-LSC-ABRM-HC-IHI-FST-A-RD	10.0	10	25	.032	18 TW

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

### Tray components include:

- Power-injectable Cook Spectrum® central venous catheter with hydrophilic coating
- Safe-T-J® double flexible-tipped wire guide with cm markings
- EchoTip® echogenic introducer needle
- Catheter insertion checklist
- BakSnap® safety syringe
- Prefilled sodium chloride syringes
- 25 gage needle
- 22 gage needle
- FEP catheter introducer needle
- Transducer tubing
- Dilator
- Lidocaine and lidocaine label
- Filter straw
- Tinted Chloraprep® One-Step
- Full-body fenestrated drape with clear window
- Gauze sponges
- Disposable safety scalpel
- Straight needle with suture or curved needle with suture and needle holder
- Disposable syringes
- MicroCLAVE® injection caps
- Locking sharps container
- Movable suture wing
- Cap, mask with face shield, and gown

**Spectrum tray configured to meet AHRQ, CDC, IHI, Joint Commission, OSHA and SHEA guidelines for patient safety.**






# Ciaglia Blue Rhino® G2

ADVANCED PERCUTANEOUS TRACHEOSTOMY INTRODUCER

Blue Rhino® G2 is a percutaneous tracheostomy introducer that enables full dilatation of the stoma in just one pass. The device's crosshatch handle helps enhance the operator's grip to improve procedural control, while longitudinal grooves on the distal surface help facilitate a less forceful insertion. With these unique design features and a comprehensive procedural tray, Blue Rhino G2 provides you with a simple, effective solution for bedside percutaneous dilatational tracheostomy (PDT) procedures.

Compared to surgical tracheostomy, bedside PDT performed with the Blue Rhino G2 offers substantial cost savings due to the elimination of high OR expenses and anesthesia fees as well as reduced physician time.

**Cost Difference of PDT vs. Surgical Tracheostomy (ST)**





Order Number	Reference Part Number	Catheter Fr
G53177	C-PTISY-100-HC-G	24.0, 26.0, 28.0
G53176	C-PTISY-100-FLEX-HC-G	21.0, 24.0, 27.0
G55629	C-PTISY-300-HC-G	21.0, 24.0, 27.0
G55630	C-PTISY-350-HC-G*	21.0, 24.0, 27.0
With VersaTube Tapered Tracheostomy Tubes		
G55658	C-PTISY-300-HC-G-VT7	21.0, 24.0, 27.0
G55659	C-PTISY-300-HC-G-VT8	21.0, 24.0, 27.0
G55660	C-PTISY-350-HC-G-VT9*	21.0, 24.0, 27.0
With Shiley Percutaneous Tracheostomy Tubes		
G53178	C-PTISY-100-HC-G-PERC6	24.0, 26.0, 28.0
G53179	C-PTISY-100-HC-G-PERC8	24.0, 26.0, 28.0

\* Does not include 15 gage, 5 cm introducer needle, and 14 Fr dilator is 6.5 cm

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

### Tray components include:

- Ciaglia Blue Rhino G2 percutaneous tracheostomy dilator with preloaded guiding catheter
- .052 inch diameter wire guide
- 15 gage, 5 cm introducer needle
- 15 gage, 7 cm introducer needle
- 15 gage, 7 cm FEP sheath needle
- 14 Fr, 4.5 cm dilator
- Tracheostomy tube loading dilator
- Disposable #15 scalpel
- Curved mosquito clamps
- 22 gage infiltration needle
- 25 gage infiltration needle
- ChlorPrep® One-Step
- Lidocaine
- Double-swivel connector
- Suture with needle
- Needle driver
- Dale® adjustable tracheostomy tube holder
- Needle holder cup
- Disposable syringes
- Fenestrated drape
- Kendall® Excilon® split tracheostomy dressing
- Lubricating jelly
- Gauze sponges

### Also available with:

- Tracheostomy tube



# Dolphin BT™

CIAGLIA BALLOON-ASSISTED TRACHEOSTOMY INTRODUCER

Dolphin BT™ is the first device to combine stoma formation and tracheal tube insertion into a simple and efficient one-step procedure: balloon-assisted tracheostomy (BAT). It was designed specifically to maximize operator control and help prevent trauma to the anterior tracheal wall and surrounding structures. Like percutaneous tracheostomy, BAT can be performed at the patient's bedside and yields comparable cost savings.

## Balloon-Assisted Tracheostomy - Procedural Benefits

- Reduced risk of trauma
- Prevention of airway leakage and blood aerosolization
- Improved preservation of tissue elasticity
- Efficiency and convenience for clinicians



Order Number	Reference Part Number	Loading Dilator Fr	Tracheostomy Tube ISO size
G36107	C-PTBSY-2400	24	
G36108	C-PTBSY-2600	26	
G36109	C-PTBSY-2700	27	
G36110	C-PTBSY-2800	28	
G49506	C-PTBSY-3000	30	
With VersaTube™ Tapered Tracheostomy Tube*			
G55594	C-PTBSY-2400-VT8	24	8
G55596	C-PTBSY-2700-VT9	27	9
* 21 Fr Dolphin BT tray with 7 mm VersaTube available in 2011.			Tracheostomy Tube Jackson size
With Shiley Percutaneous Cuffed Tracheostomy Tube			
G36118	C-PTBSY-2600-PERC6	26	6
G49507	C-PTBSY-2800-PERC8	28	8

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

### Tray components include:

- Dolphin tracheostomy introducer
- Cook Inflation Device
- 18 gage introducer EchoTip® needle
- 18 gage FEP sheath needle
- Wire guide with positioning marks
- 14 Fr dilator
- Disposable #15 safety scalpel
- Curved mosquito clamps
- 22 gage infiltration needle
- 25 gage infiltration needle
- Chlorhexidine/alcohol prep solution
- Lidocaine epinephrine 1.5%

- Double-swivel connector
- Suture with needle
- Needle holder cup
- Disposable syringe
- Full-body drape with clear plastic window (not pictured)
- Lubricating jelly (not pictured)
- Gauze sponges (not pictured)

### Also available with:

- Tracheostomy tube




## Tiger 2™

SELF-ADVANCING NASAL JEJUNAL FEEDING TUBE

Tiger 2™ is an innovative self-advancing nasal jejunal feeding tube that allows peristalsis to pull it directly into the jejunum, providing quick enteral access for the delivery of nutrition and medication. It can be placed by nurses at the bedside, avoiding costly endoscopy and fluoroscopy procedures.

In a recent clinical trial,<sup>10</sup> Tiger 2 achieved a postpyloric **placement success rate of 90%**, and the average clinician time required to place each tube was less than 10 minutes.


Cilia-like flaps allow the feeding tube to self advance, aid in reducing the incidence of dislodgment and eliminate the need to reposition the patient.



Distance markings every 10 cm from 40-100 cm provide visual confirmation of tube position (radiopaque for enhanced x-ray visualization).

Order Number	Reference Part Number	Fr	Length cm	Sideports
<b>Tiger 2</b>				
G50415	C-NJFT-14-U	14.0	155	5
		Diameter inch	Length cm	Tip Configuration
<b>Torque Cable (optional)</b>				
G30124	TSBD-65-170-NJFT	.065	170	straight

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.



Five large sideports help to prevent the tube from clogging.

14 Fr diameter optimizes flow of medications and thicker, fiber-containing formulas.



# DIFFICULT AIRWAY PRODUCTS

## Intubation

Cook offers a variety of airway catheters for blind or bronchoscopic-guided intubation. All intubation catheters feature blunt tips that are atraumatic to internal structures.

## Exchange

Unlike traditional stylets used to facilitate tube exchange, Cook’s airway exchange catheters are hollow to permit oxygen delivery if necessary. They are also available in different size options and materials to accommodate varying clinical needs.

## Controlled Extubation


Even with well-designed protocols, postextubation airway failure remains a significant risk in all patient populations. Cook’s airway exchange catheter can be used to maintain continuous airway access for oxygen delivery and facilitate a successful reintubation if required.

## Emergency Airway

When traditional treatment methods have failed, Cook provides life-saving solutions with emergency cricothyrotomy and transtracheal catheters to quickly and safely restore patients’ airways.



# AIRWAY EXCHANGE CATHETERS




### Components

- Radiopaque catheter
- Rapi-Fit® Adapter (15 mm connector)
- Rapi-Fit Adapter (Luer lock connector)

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter ID mm	For Placement of ET Tubes* with ID mm
For Use with Single-Lumen Endotracheal Tubes					
G07833	C-CAE-8.0-45	8.0	45	1.6	3 or larger
G06732	C-CAE-11.0-83	11.0	83	2.3	4 or larger
G07873	C-CAE-14.0-83	14.0	83	3.0	5 or larger
G05880	C-CAE-19.0-83	19.0	83	3.4	7 or larger

\*Endotracheal tube not included  
Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.



**Extra Firm with Soft Tip  
Components**

- Radiopaque catheter
- Rapi-Fit® Adapter (15 mm connector)
- Rapi-Fit Adapter (Luer lock connector)

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter ID mm	For Placement of ET Tubes* with ID mm
For Use with Double-Lumen Endotracheal Tubes					
G36401	C-CAE-11.0-100-DLT-EF-ST	11.0	100	2.3	4 or larger
G36402	C-CAE-14.0-100-DLT-EF-ST	14.0	100	3.0	5 or larger


\*Endotracheal tube not included.  
Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.





# Frova

## INTUBATING INTRODUCER



### Components

- Radiopaque catheter
- Rapi-Fit® Adapter (15 mm connector)
- Rapi-Fit Adapter (Luer lock connector)

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter ID mm	For Placement of ET Tubes* with ID mm
With Stiffening Cannula					
G12591	C-CAE-8.0-35-FII	8.0	35	1.6	3 or larger
G11915	C-CAE-14.0-65-FII	14.0	65	3	6 or larger
Without Stiffening Cannula					
G13307	C-CAE-14.0-65-FI	14.0	65	3	6 or larger
Without Stiffening Cannula and Without Rapi-Fit Adapter**					
G27282	C-CAE-14.0-65-FIC	14.0	65	3	6 or larger

\*Endotracheal tube not included.


\*\*Supplied as 10 catheters only.

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.



# Aintree

## INTUBATION CATHETER



### Components

- Radiopaque catheter
- Rapi-Fit® Adapter (15 mm connector)
- Rapi-Fit Adapter (Luer lock connector)

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter ID mm	For Placement of ET Tubes* with ID mm
For Use with Laryngeal Mask Airways/Supraglottic Airway Devices					
G10789	C-CAE-19.0-56-AIC	19.0	56	4.7	7 or larger

\*Endotracheal tube not included.

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.



# Enk

## OXYGEN FLOW MODULATOR SET



### Components

- Catheter needle
- Syringe
- Enk Oxygen Flow Modulator
- Emergency Transtracheal Airway Catheter tubing

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter ID mm	Needle gage
<b>Sets</b>					
G12398	C-EFMS-100	6.0	7.5	2	15
G13145	C-EFMS-101	6.0	5.0	2	15
<b>Emergency Transtracheal Airway Catheters</b>					
G10338	C-DTJV-6.0-5.0-BTT	6.0	5.0	2	15
G09413	C-DTJV-6.0-7.5-BTT	6.0	7.5	2	15

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.



# Melker

## EMERGENCY CRICOTHYROTOMY CATHETER



### Universal Set - Seldinger and Surgical Technique

#### Components

- Introducer needle
- FEP catheter
- Syringe
- Scalpel
- Wire guide
- Curved dilator
- Airway catheter
- #11 safety scalpel
- Trousseau dilator
- Tracheal hook
- Blunt, curved dilator
- Trach tape
- Tinted Chloraprep® One-Step
- Gauze sponges
- Needle cup holder

Order Number	Reference Part Number	Catheter Length cm	Catheter ID mm
G32193	C-TCCSB-500-UNI	9	5

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.



## CENTESIS & DRAINAGE

Cook's full line of catheters for bedside drainage procedures includes a variety of sizes and design configurations. Some devices can be used to drain both air and fluid, minimizing the number of thoracostomy tools your hospital must keep in stock. To accommodate physician preference, options are available for both the trocar and Seldinger insertion technique. Additionally, Cook catheters are packaged in comprehensive procedural trays with all the necessary components for successful and efficient bedside placement.



# Thal-Quick

## CHEST TUBE INTRODUCER





Single Lumen



Double Lumen

Order Number	Reference Part Number	Fr	Length cm	Sideports
<b>Single Lumen Trays</b>				
G06791	C-TQTSY-800	8.0	18	3
G06885	C-TQTSY-1000	10.0	20	3
G05464	C-TQTSY-1200	12.0	22	3
G07090	C-TQTSY-1400	14.0	41	4
G05462	C-TQTSY-1600	16.0	41	4
G05461	C-TQTSY-1800	18.0	41	4
G05463	C-TQTSY-2000	20.0	41	4
G04220	C-TQTSY-2400	24.0	41	4
G04226	C-TQTSY-2800	28.0	41	4
G07242	C-TQTSY-3200	32.0	41	4
G07244	C-TQTSY-3600	36.0	41	4
<b>Double Lumen Tray</b>				
G36225	C-DLCTSY-2400	24.0	40	4 large, 4 small

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

## Thal-Quick Chest Tube Trays

### Components

- Scalpel
- EchoTip® echogenic introducer needle
- Amplatz Extra-Stiff Wire Guide
- Dilators
- Chest tube inserter
- Thal-Quick Chest Tube
- 25 gage needle
- 22 gage needle
- Lidocaine
- Chlorhexidine
- Fenestrated drape
- Gauze sponges
- Disposable syringes
- Suture with needle
- Povidone-iodine ointment
- CSR wrap
- Prep tray
- Needle holder cup



# Wayne

## PNEUMOTHORAX CATHETER



### Seldinger Set Components

- EchoTip® echogenic introducer needle
- 9 cm introducer needles
- Wire guide
- Dilator
- Pigtail catheter
- Cook Chest Drain Valve
- Connecting tube
- #11 thumb scalpel

### Trocar Set Components

- Needle obturator
- Pigtail catheter
- Cook Chest Drain Valve
- Connecting tube
- #11 thumb scalpel

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Needle gage	Needle Length cm
<b>Seldinger Set*</b>						
G11182	C-UTPT-1400-WAYNE-112497	14.0	29	19	18	5, 9
<b>Trocar Sets*</b>						
G09922	C-UTPT-1020-WAYNE	10.2	29	20	19	16.5
G07910	C-UTPT-1400-WAYNE	14.0	29	20	19	16.5

\*Sets also include three-way stopcock.

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.



# Lock

## PERICARDIOCENTESIS CATHETER SET



### Components

- Introducer needles
- Wire guide
- Dilator
- Radiopaque pigtail catheter
- Multipurpose tubing adapter

Order Number	Reference Part Number	Catheter Fr/Length cm	Catheter Sideports	Wire Guide Diameter/Length cm/inch	Wire Guide Tip Configuration	Needle gage/Length cm
G02433	C-PCS-830-LOCK	8.3/40	6	.038/70	3 mm J	18/7, 12

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.



# Fuhrman

## PLEURAL/PNEUMOPERICARDIAL DRAINAGE SET



### Components

- Introducer needle
- Wire guide
- Dilator
- Radiopaque pigtail catheter
- Three-way stopcock
- Multipurpose tubing adapter

Order Number	Reference Part Number	Catheter Fr/Length cm	Catheter Sideports	Wire Guide Diameter/Length inch/cm	Wire Guide Tip Configuration	Needle gage/Length cm
G55723	C-PPD-850-IMH	8.5/15	6	.028/40	2 mm J	18/7

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.



## PERITONEAL LAVAGE SET



### Components

- Introducer needle
- Straight safety wire guide
- Radiopaque catheter
- 25 gage needle
- 22 gage needle
- Lidocaine
- Chlorhexidine
- Fenestrated drape
- Gauze sponges
- Disposable syringes
- Disposable scalpel
- Suture with needle
- Povidone-iodine ointment
- CSR wrap
- Prep tray
- Needle holder cup

Order Number	Reference Part Number	Catheter Fr/Length cm	Catheter Sideports	Wire Guide Diameter/Length inch/cm	Needle gage/Length cm
G05268	C-PLSY-901	9.0/20	90	.038/40	18/7

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

## References

1. Waydhas C. Equipment review: intrahospital transport of critically ill patients. *Crit Care*. 1999;3(5):R83-R89.
2. Shippert RD. A study of time-dependent operating room fees and how to save \$100,000 by using time-saving products. *Am J Cosmet Surg*. 2005;22(1):25-34.
3. Schuerer DJ, Zach JE, Thomas J, et al. Effect of chlorhexidine/silver sulfadiazine-impregnated central venous catheters in an intensive care unit with a low blood stream infection rate after implementation of an educational program: a before-after trial. *Surg Infect (Larchmt)*. 2007;8(4):445-454.
4. Schuerer DJ, Mazuski JE, Buchman TJ, et al. Catheter-related bloodstream infection rates in minocycline/rifampin vs. chlorhexidine/silver sulfadiazine-impregnated central venous catheters-results of a 46 month study. *Crit Care Med*. 2008;36(12)(suppl):A199-A208. Abstract 454.
5. Warren DK, Quadir WW, Hollenbeak CS, et al. Attributable cost of catheter-associated bloodstream infections among intensive care patients in a nonteaching hospital. *Crit Care Med*. 2006;34(8):2084-2089.
6. O'Grady NP, Alexander M, Dellinger EP, et al.; Centers for Disease Control and Infection. Guidelines for the prevention of intravascular catheter-related infections. *MMWR Recomm Rep*. 2002;51(RR):1-26.
7. Bowen CP, Whitney LR, Truitt JD, et al. Comparison of safety and cost of percutaneous versus surgical tracheostomy. *Am Surg*. 2001;67(1):54-60.
8. Cobean R, Beals M, Moss C, et al. Percutaneous dilatational tracheostomy. A safe, cost-effective bedside procedure. *Arch Surg*. 1996;131(3):265-271.
9. Freeman BD, Isabella K, Cobb JP, et al. A prospective, randomized study comparing percutaneous with surgical tracheostomy in critically ill patients. *Crit Care Med*. 2001;29(5):926-930.
10. Samis AJW, Heyland D, Drover J. Evaluation of two novel strategies for post-pyloric placement of enteral feeding tubes. Poster presented at: Toronto Critical Care Medicine Symposium; October 30-November 1, 2003; Toronto, Canada.

*BakSnap is a registered trademark of DuoProSS Meditech Corporation.*

*ChlorPrep is a registered trademark of CareFusion 2200.*


*MicroCLAVE is a registered trademark of ICU Medical, Inc.*

*Dale is a registered trademark of Dale Medical Products, Inc.*

*Kendall and Excilon are registered trademarks of Tyco Healthcare Group LP.*







AORTIC INTERVENTION  
CRITICAL CARE  
ENDOSCOPY  
INTERVENTIONAL RADIOLOGY  
LEAD MANAGEMENT  
PERIPHERAL INTERVENTION  
SURGERY  
UROLOGY  
WOMEN'S HEALTH



#### Customer Service Centers

**EMEA:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
**Distributors:** +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)  
**Austria:** +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)  
**Belgium:** +32 27001633, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)  
**Denmark:** +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)  
**France:** +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
**Germany:** +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
**Hungary:** +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)  
**Ireland:** +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)  
**Italy:** +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
**Netherlands:** +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)  
**Norway:** +47 23162968, [no.orders@cookmedical.com](mailto:no.orders@cookmedical.com)  
**Poland:** +48 223060159, [pl.orders@cookmedical.com](mailto:pl.orders@cookmedical.com)  
**Spain:** +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)  
**Sweden:** +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)  
**Switzerland - French:** +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
**Switzerland - Italian:** +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
**Switzerland - German:** +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
**United Kingdom:** +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)  
**[www.cookmedical.com](http://www.cookmedical.com)**

**Americas:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335,  
E-mail: [orders@cookmedical.com](mailto:orders@cookmedical.com)

**Australia:**  
Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283,  
E-mail: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)

