

# Band with confidence.



## 4,6,10 Shooter®

SAEED® MULTI-BAND LIGATOR

After leading the market for more than a decade, the reliable Saeed Multi-Band Ligator is now available with Universal-fit barrels that streamline your order information and simplify inventory control. And, by using bands that demonstrate lower slippage rates when compared to other bands on the market<sup>1</sup>, the Saeed 4, 6, 10 Shooter can potentially reduce the risk of rebleeding. We've taken reliability and streamlined it.

<sup>1</sup> "In-vitro bench studies have shown differences in the resting inner diameters and radial retraction force of bands on different ligation devices. The larger resting inner diameter and a significantly lower retraction force in neoprene bands, compared to rubber bands predisposed these bands to early slippage." J. M. Polski, E. M. Brunt, Z. A. Saeed (2001). Chronology of Histological Changes after Band Ligation of Esophageal Varices in Humans. *Endoscopy*, 33, 443-447.

<sup>2</sup> Ibid, p.446



[www.cookmedical.com](http://www.cookmedical.com)



## 4,6,10 Shooter\*

SAEED® MULTI-BAND LIGATOR

Used to endoscopically ligate esophageal varices at or above the gastroesophageal junction or to ligate internal hemorrhoids.

Supplied non-sterile and is disposable - intended for single use only.

- Universal-fit barrels streamline ligation procedures and inventory control
- Opti-Vu® Barrel features two-string trigger cord and recessed bands – providing maximum endoscopic view while maintaining suction
- Clear band signals when only one band remains to be deployed
- Wide array of barrel sizes to meet a full range of procedural needs

\*The larger resting inner diameter and a significantly lower retraction force in neoprene bands, compared to rubber bands predisposed these bands to early slippage.<sup>2</sup>

Order Number	Reference Part Number	Endoscope Outer Diameter mm	Number of Bands	Trigger Cord Length cm	Order Number	Reference Part Number	Endoscope Outer Diameter mm	Number of Bands	Trigger Cord Length cm
G31915	MBL-U-4	8.6 - 11.3	4	122	G25423	MBL-4-XS	8.6 - 9.2	4	122
G31917	MBL-U-6	8.6 - 11.3	6	122	G25424	MBL-6-XS	8.6 - 9.2	6	122
G31919	MBL-U-10	8.6 - 11.3	10	122	G25504	MBL-10-XS	8.6 - 9.2	10	122
G31916	MBL-U-4-F	8.6 - 11.3	4	142	G25122	MBL-4-XL	11 - 14	4	122
G31918	MBL-U-6-F	8.6 - 11.3	6	142	G25119	MBL-6-XL	11 - 14	6	122
G31920	MBL-U-10-F	8.6 - 11.3	10	142	G25220	MBL-4-XL-C	11 - 14	4	190
G24654	MBL-6-OV	9.5 - 11.5	6	122	G25221	MBL-6-XL-C	11 - 14	6	190
G25503	MBL-6-OV-XS	8.6 - 9.2	6	122	G22690	MBL-4-F	9.5 - 13	4	142
G22613	MBL-4	9.5 - 13	4	122	G22692	MBL-6-F	9.5 - 13	6	142
G22553	MBL-6	9.5 - 13	6	122	G22691	MBL-10-F	9.5 - 11.5	10	142
G22614	MBL-10	9.5 - 11.5	10	122	G24655	MBL-6-F-OV	9.5 - 13	6	142

Each kit includes Opti-Vu barrel w/ ligation bands and trigger cord, loading catheter, precision control handle and irrigation adapter. Minimum accessory channel 2.8 mm.

**Caution:** This product contains natural rubber latex, which may cause an allergic reaction.

**Note:** The new universal-fit barrels are not compatible with XL scopes.

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

### Customer Service

**EMEA:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
 Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)  
 Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)  
 Belgium: +32 27001633, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)  
 Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)  
 Finland: +358 972519996, [fi.orders@cookmedical.com](mailto:fi.orders@cookmedical.com)  
 France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)  
 Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)  
 Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)  
 Norway: +47 23162968, [no.orders@cookmedical.com](mailto:no.orders@cookmedical.com)  
 Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)  
 Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)  
 Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

[www.cookmedical.com](http://www.cookmedical.com)

**Americas:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335  
 E-mail: [orders@cookmedical.com](mailto:orders@cookmedical.com)

**Australia:**  
 Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283  
 E-mail: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)

