

Cross lesions like never before.



Advance[®] 18LP

LOW-PROFILE PTA BALLOON DILATATION CATHETER

Broad offering

Advance 18LP is available in nearly 200 size, length, and shaft configuration options.

Trackability

A flexible tip enhances your ability to navigate tortuous vasculature.

Fast Rewrap

Advanced rewrap technology provides fast balloon deflation so you can quickly restore blood flow.



The following products are compatible with Advance 18LP.

ACCESS




Micropuncture®
Introducer Set

TARGET




Flexor® Ansel Guiding Sheath
with High-Flex Dilator




Slip-Cath®
Angiographic Catheters



Approach® CTO
Microwire Guide



CXI®
Support Catheter



Approach® Hydro ST
Microwire Guide

TREAT



Zilver 518®
Vascular Self-Expanding Stent


ADVANCE 18LP LOW-PROFILE PTA DILATATION CATHETER

Order Number	Reference Part Number	Compatible Sheath Fr	Nominal Pressure atm	Rated Burst Pressure atm
80 cm Catheter Length				
Over the wire				
G30892	PTA4-18-80-4-1.5	4	8	14
G50225	PTA4-18-80-4-2	4	8	14
G30893	PTA4-18-80-4-2.5	4	8	14
G50227	PTA4-18-80-4-4	4	8	14
G50229	PTA4-18-80-4-6	4	8	14
G50231	PTA4-18-80-4-8	4	8	14
G50232	PTA4-18-80-4-10	4	8	14
G30894	PTA4-18-80-4-12	4	8	12
G30895	PTA4-18-80-4-15	4	8	12
G30896	PTA4-18-80-4-17	4	8	12
G30897	PTA4-18-80-4-20	4	8	12
G30898	PTA4-18-80-5-1.5	4	8	14
G30899	PTA4-18-80-5-2.5	4	8	14
G50234	PTA4-18-80-5-2	4	8	14
G50236	PTA4-18-80-5-4	4	8	14
G50238	PTA4-18-80-5-6	4	8	14
G50240	PTA4-18-80-5-8	4	8	14
G50242	PTA4-18-80-5-10	4	8	14
G30900	PTA4-18-80-5-12	4	8	12
G30901	PTA4-18-80-5-15	4	8	12
G30902	PTA4-18-80-5-17	4	8	12
G30903	PTA4-18-80-5-20	4	8	12
G30904	PTA4-18-80-6-1.5	4	8	14
G50244	PTA4-18-80-6-2	4	8	14
G30905	PTA4-18-80-6-2.5	4	8	14
G50246	PTA4-18-80-6-4	4	8	14
G50248	PTA4-18-80-6-6	4	8	14
G50250	PTA4-18-80-6-8	4	8	14
G50252	PTA4-18-80-6-10	4	8	14
G30906	PTA4-18-80-6-12	5	8	12
G30907	PTA4-18-80-6-15	5	8	12
G30908	PTA4-18-80-6-17	5	8	12
G30909	PTA4-18-80-6-20	5	8	12
G30910	PTA4-18-80-7-1.5	4	8	14
G50254	PTA4-18-80-7-2	4	8	14
G30911	PTA4-18-80-7-2.5	4	8	14
G50256	PTA4-18-80-7-4	4	8	14
G50258	PTA4-18-80-7-6	4	8	14
G50260	PTA4-18-80-7-8	4	8	14
G50263	PTA4-18-80-7-10	4	8	14
G30912	PTA4-18-80-7-12	5	8	10
G30913	PTA4-18-80-7-15	5	8	10
G30914	PTA4-18-80-7-17	5	8	10

Order Number	Reference Part Number	Compatible Sheath Fr	Nominal Pressure atm	Rated Burst Pressure atm
80 cm Catheter Length				
Over the wire				
G30915	PTA4-18-80-7-20	5	8	10
G30916	PTA4-18-80-8-1.5	5	8	14
G50265	PTA4-18-80-8-2	5	8	14
G30917	PTA4-18-80-8-2.5	5	8	14
G50267	PTA4-18-80-8-4	5	8	14
G50269	PTA4-18-80-8-6	5	8	14
G50271	PTA4-18-80-8-8	5	8	14
G50273	PTA4-18-80-8-10	5	8	14
G30918	PTA4-18-80-8-12	6	8	10
G30919	PTA4-18-80-8-15	6	8	10
G30920	PTA4-18-80-8-17	6	8	10
G30921	PTA4-18-80-8-20	6	8	10
G30922	PTA4-18-80-9-2	5	8	10
G30923	PTA4-18-80-9-4	5	8	10
G30924	PTA4-18-80-9-6	5	8	10
G30925	PTA4-18-80-9-8	5	8	10
G30926	PTA4-18-80-9-10	6	8	10
G30927	PTA4-18-80-10-2	6	8	10
G30928	PTA4-18-80-10-4	6	8	10
G30929	PTA4-18-80-10-6	6	8	10
G30930	PTA4-18-80-10-8	7	8	10
G30931	PTA4-18-80-10-10	7	8	10

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Reference Part Number Key



Order Number	Reference Part Number	Compatible Sheath Fr	Nominal Pressure atm	Rated Burst Pressure atm
150 cm Catheter Length				
Over the wire				
G30932	PTA4-18-150-2-2	4	8	14
G30933	PTA4-18-150-2-4	4	8	14
G30934	PTA4-18-150-2-6	4	8	14
G30935	PTA4-18-150-2-8	4	8	14
G34166	PTA4-18-150-2-10	4	8	14
G30937	PTA4-18-150-2-12	4	8	14
G30938	PTA4-18-150-2-15	4	8	14
G30939	PTA4-18-150-2-17	4	8	14
G30940	PTA4-18-150-2-20	4	8	14
G30941	PTA4-18-150-3-1.5	4	8	14
G34318	PTA4-18-150-3-2	4	8	14
G30942	PTA4-18-150-3-2.5	4	8	14
G34319	PTA4-18-150-3-4	4	8	14
G34320	PTA4-18-150-3-6	4	8	14
G34321	PTA4-18-150-3-8	4	8	14
G34322	PTA4-18-150-3-10	4	8	14
G30943	PTA4-18-150-3-12	4	8	12
G30944	PTA4-18-150-3-15	4	8	12
G30945	PTA4-18-150-3-17	4	8	12
G30946	PTA4-18-150-3-20	4	8	12
G30947	PTA4-18-150-4-1.5	4	8	14
G34323	PTA4-18-150-4-2	4	8	14
G30948	PTA4-18-150-4-2.5	4	8	14
G34324	PTA4-18-150-4-4	4	8	14
G34325	PTA4-18-150-4-6	4	8	14
G34326	PTA4-18-150-4-8	4	8	14
G34327	PTA4-18-150-4-10	4	8	14
G30949	PTA4-18-150-4-12	4	8	12
G30950	PTA4-18-150-4-15	4	8	12
G30951	PTA4-18-150-4-17	4	8	12
G30952	PTA4-18-150-4-20	4	8	12
G30953	PTA4-18-150-5-1.5	4	8	14
G34328	PTA4-18-150-5-2	4	8	14
G30954	PTA4-18-150-5-2.5	4	8	14
G34329	PTA4-18-150-5-4	4	8	14
G34330	PTA4-18-150-5-6	4	8	14
G34331	PTA4-18-150-5-8	4	8	14
G34332	PTA4-18-150-5-10	4	8	14
G30955	PTA4-18-150-5-12	4	8	12
G30956	PTA4-18-150-5-15	4	8	12
G34167	PTA4-18-150-5-17	4	8	12
G30958	PTA4-18-150-5-20	4	8	12
G30959	PTA4-18-150-6-1.5	4	8	14

Order Number	Reference Part Number	Compatible Sheath Fr	Nominal Pressure atm	Rated Burst Pressure atm
150 cm Catheter Length				
Over the wire				
G34333	PTA4-18-150-6-2	4	8	14
G30960	PTA4-18-150-6-2.5	4	8	14
G34334	PTA4-18-150-6-4	4	8	14
G34335	PTA4-18-150-6-6	4	8	14
G34336	PTA4-18-150-6-8	4	8	14
G34337	PTA4-18-150-6-10	4	8	14
G30961	PTA4-18-150-6-12	5	8	12
G30962	PTA4-18-150-6-15	5	8	12
G30963	PTA4-18-150-6-17	5	8	12
G30964	PTA4-18-150-6-20	5	8	12
G30965	PTA4-18-150-7-1.5	4	8	14
G34338	PTA4-18-150-7-2	4	8	14
G30966	PTA4-18-150-7-2.5	4	8	14
G34339	PTA4-18-150-7-4	4	8	14
G34340	PTA4-18-150-7-6	4	8	14
G34341	PTA4-18-150-7-8	4	8	14
G34342	PTA4-18-150-7-10	4	8	14
G30967	PTA4-18-150-7-12	5	8	10
G30968	PTA4-18-150-7-15	5	8	10
G30969	PTA4-18-150-7-17	5	8	10
G30970	PTA4-18-150-7-20	5	8	10
G30971	PTA4-18-150-8-1.5	5	8	14
G34343	PTA4-18-150-8-2	5	8	14
G30972	PTA4-18-150-8-2.5	5	8	14
G34344	PTA4-18-150-8-4	5	8	14
G34345	PTA4-18-150-8-6	5	8	14
G34346	PTA4-18-150-8-8	5	8	14
G34347	PTA4-18-150-8-10	5	8	14
G30973	PTA4-18-150-8-12	6	8	10
G30974	PTA4-18-150-8-15	6	8	10
G30975	PTA4-18-150-8-17	6	8	10
G30976	PTA4-18-150-8-20	6	8	10
G30977	PTA4-18-150-9-2	5	8	10
G30978	PTA4-18-150-9-4	5	8	10
G30979	PTA4-18-150-9-6	5	8	10
G30980	PTA4-18-150-9-8	5	8	10
G30981	PTA4-18-150-9-10	6	8	10

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Order Number	Reference Part Number	Compatible Sheath Fr	Nominal Pressure atm	Rated Burst Pressure atm
170 cm Catheter Length				
Peripheral exchange				
G30815	PTAX4-18-170-2-2	5	8	14
G30816	PTAX4-18-170-2-4	5	8	14
G30817	PTAX4-18-170-2-6	5	8	14
G30818	PTAX4-18-170-2-8	5	8	14
G30819	PTAX4-18-170-2-10	5	8	14
G30820	PTAX4-18-170-2-12	5	8	14
G30821	PTAX4-18-170-2-15	5	8	14
G30822	PTAX4-18-170-2-17	5	8	14
G30823	PTAX4-18-170-2-20	5	8	14
G30824	PTAX4-18-170-3-2	5	8	14
G30825	PTAX4-18-170-3-4	5	8	14
G30826	PTAX4-18-170-3-6	5	8	14
G30827	PTAX4-18-170-3-8	5	8	14
G30828	PTAX4-18-170-3-10	5	8	14
G30829	PTAX4-18-170-3-12	5	8	12
G30830	PTAX4-18-170-3-15	5	8	12
G30831	PTAX4-18-170-3-17	5	8	12
G30832	PTAX4-18-170-3-20	5	8	12
G30833	PTAX4-18-170-4-2	5	8	14
G30834	PTAX4-18-170-4-4	5	8	14
G30835	PTAX4-18-170-4-6	5	8	14
G30836	PTAX4-18-170-4-8	5	8	14
G30837	PTAX4-18-170-4-10	5	8	14
G30838	PTAX4-18-170-4-12	5	8	12
G30839	PTAX4-18-170-4-15	5	8	12
G30840	PTAX4-18-170-4-17	5	8	12
G30841	PTAX4-18-170-4-20	5	8	12
G30842	PTAX4-18-170-5-2	5	8	14
G30843	PTAX4-18-170-5-4	5	8	14
G30844	PTAX4-18-170-5-6	5	8	14
G30845	PTAX4-18-170-5-8	5	8	14
G30846	PTAX4-18-170-5-10	5	8	14
G30847	PTAX4-18-170-5-12	5	8	12
G30848	PTAX4-18-170-5-15	5	8	12
G30849	PTAX4-18-170-5-17	5	8	12
G30850	PTAX4-18-170-5-20	5	8	12
G30851	PTAX4-18-170-6-2	5	8	14

Order Number	Reference Part Number	Compatible Sheath Fr	Nominal Pressure atm	Rated Burst Pressure atm
170 cm Catheter Length				
Peripheral exchange				
G30852	PTAX4-18-170-6-4	5	8	14
G30853	PTAX4-18-170-6-6	5	8	14
G30854	PTAX4-18-170-6-8	5	8	14
G30855	PTAX4-18-170-6-10	5	8	14
G30856	PTAX4-18-170-6-12	5	8	12
G30857	PTAX4-18-170-6-15	5	8	12
G30858	PTAX4-18-170-6-17	5	8	12
G30859	PTAX4-18-170-6-20	5	8	12
G30860	PTAX4-18-170-7-2	5	8	14
G30861	PTAX4-18-170-7-4	5	8	14
G30862	PTAX4-18-170-7-6	5	8	14
G30863	PTAX4-18-170-7-8	5	8	14
G30890	PTAX4-18-170-7-10	5	8	14
G30865	PTAX4-18-170-7-12	5	8	10
G30866	PTAX4-18-170-7-15	5	8	10
G30867	PTAX4-18-170-7-17	5	8	10
G30868	PTAX4-18-170-7-20	5	8	10

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001633, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 972519996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
E-mail: orders@cookmedical.com

Australia:
Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283
E-mail: cau.custserv@cookmedical.com

