


QUICK REFERENCE GUIDE

SURGICAL MINI BIOPSY FORCEPS




Cup Diameter
1.8 mm

Cup Length
2.46 mm

Open Cup Width
5.91 mm

SURGICAL BIOPSY FORCEPS




Cup Diameter
2.4 mm

Cup Length
2.95 mm

Open Cup Width
7.26 mm

SERRATED BIOPSY FORCEPS



Cup Diameter
2.4 mm

Cup Length
2.99 mm

Open Cup Width
7.58 mm

SERRATED LARGE BIOPSY FORCEPS



Cup Diameter
2.4 mm

Cup Length
3.68 mm

Open Cup Width
8.95 mm

SERRATED JUMBO BIOPSY FORCEPS




Cup Diameter
2.8 mm

Cup Length
3.93 mm

Open Cup Width
8.98 mm

SERRATED MAX BIOPSY FORCEPS



Cup Diameter
3.3 mm

Cup Length
3.98 mm

Open Cup Width
9.98 mm

HOT BIOPSY FORCEPS



Cup Diameter
2.4 mm

Cup Length
2.99 mm

Open Cup Width
7.58 mm

HOT SERRATED BIOPSY FORCEPS



Cup Diameter
2.4 mm

Cup Length
2.99 mm

Open Cup Width
7.58 mm



Captura[®]

BIOPSY FORCEPS

QUICK REFERENCE ORDER INFORMATION

Not all Reference Part Numbers are available in all markets.

Reference Part Number	Order Number	Reference Part Number	Order Number
DBF-1.8-S	G47742	DBF-2.4SL-230SP-S	G53641
DBF-1.8-160-S	G53006	DBF-2.4SL-230SP-20-S	G56087
DBF-2.4-160-S	G31577	DBF-2.4SL-230SP-40-S	G56076
DBF-2.4-160-20-S	G56061	DBF-2.8SJ-160-S	G53628
DBF-2.4-160SP-S	G31580	DBF-2.8SJ-160-20-S	G56081
DBF-2.4-160SP-20-S	G56071	DBF-2.8SJ-160-40-S	G56083
DBF-2.4-230-S	G31576	DBF-2.8SJ-160SP-S	G53650
DBF-2.4-230-20-S	G56065	DBF-2.8SJ-160SP-20-S	G56062
DBF-2.4-230SP-S	G31585	DBF-2.8SJ-160SP-40-S	G56064
DBF-2.4-230SP-20-S	G56074	DBF-2.8SJ-230-S	G53642
DBF-2.4SN-160-S	G56122	DBF-2.8SJ-230-20-S	G56085
DBF-2.4SN-160-20-S	G56123	DBF-2.8SJ-230-40-S	G56086
DBF-2.4SN-160-40-S	G56124	DBF-2.8SJ-230SP-S	G53042
DBF-2.4SN-160SP-S	G56128	DBF-2.8SJ-230SP-20-S	G56066
DBF-2.4SN-160SP-20-S	G56129	DBF-2.8SJ-230SP-40-S	G56067
DBF-2.4SN-160SP-40-S	G56130	DBF-2.8SJN-230-20-S	G56721
DBF-2.4SN-230-S	G56125	DBF-2.8SJN-230-40-S	G56722
DBF-2.4SN-230-20-S	G56126	DBF-2.8SJN-230SP-20-S	G56723
DBF-2.4SN-230-40-S	G56127	DBF-2.8SJN-230SP-40-S	G56724
DBF-2.4SN-230SP-S	G56131	DBF-3.3SM-160-S	G54166
DBF-2.4SN-230SP-20-S	G56132	DBF-3.3SM-160-20-S	G56054
DBF-2.4SN-230SP-40-S	G56133	DBF-3.3SM-160SP-S	G54168
DBF-2.4SL-160-S	G53630	DBF-3.3SM-160SP-20-S	G56053
DBF-2.4SL-160-20-S	G56057	DBF-3.3SM-230-S	G54167
DBF-2.4SL-160-40-S	G56056	DBF-3.3SM-230-20-S	G56058
DBF-2.4SL-160SP-S	G53638	DBF-3.3SM-230SP-S	G54170
DBF-2.4SL-160SP-20-S	G56070	DBF-3.3SM-230SP-20-S	G56055
DBF-2.4SL-160SP-40-S	G56073	DBF-3.3SM-230SP-40-S	G56060
DBF-2.4SL-230-S	G53633	HDBF-2.4-230-S	G31583
DBF-2.4SL-230-20-S	G56052	HDBF-2.4SN-230-S	G56135
DBF-2.4SL-230-40-S	G56059		



www.cookmedical.com

COOK ENDOSCOPY
 4900 Bethania Station Road, Winston-Salem, NC 27105 U.S.A.
 Phone: 336 744-0157, Toll Free: (USA) 800 457-4500, Fax: 336 744-5231

COOK IRELAND LTD.
 O'Halloran Road, National Technology Park, Limerick, IRELAND
 Phone: +353 61 334440, Fax: +353 61 334441

© COOK 2011 L-12712/0211

AORTIC INTERVENTION

CRITICAL CARE

ENDOSCOPY

INTERVENTIONAL RADIOLOGY

LEAD MANAGEMENT

PERIPHERAL INTERVENTION

SURGERY

UROLOGY

WOMEN'S HEALTH