


COOK[®]
MEDICAL

Celect[®]

VENA CAVA FILTER



Safe, effective, dependable pulmonary embolism prevention with optional retrieval.

LEAVE IT OR RETRIEVE IT

The Celect functions well as both a permanent vena cava filter and a retrievable device.

ACCURATE PLACEMENT

Secondary strut design centers the filter with minimal ingrowth.

SECURE FIXATION

Specially designed anchors for secure, atraumatic caval fixation.



Now available with NavAlign™ delivery.

- Hemostasis valve
- Added safety mechanism with one-push filter delivery
- Femoral and jugular delivery options
- Femoral option includes tactile bump indicating filter is fully advanced
- NavAlign multipurpose dilator includes radiopaque sizing bands and flushing sideports

Global Product Number	Order Number	Filter Maximum Diameter mm	Filter Length mm	Introducer Sheath Fr	Introducer Sheath Length cm
Femoral					
G52914	IGTCFS-65-1-FEM-CELECT	30	48	7.0	65
Jugular					
G52909	IGTCFS-65-1-JUG-CELECT	30	48	7.0	65
UniSet					
G52915	IGTCFS-65-1-UNI-CELECT	30	48	7.0	65

Global Product Number	Order Number	Retrieval Loop System Catheter Fr	Retrieval Loop System Length cm	Retrieval Sheath Fr	Retrieval Sheath Length cm
Vena Cava Filter Retrieval Set (Filter Not Included)*					
G13287	GTRS-200-RB	6.3	80	11.0	60

*Set contains the following: coaxial sheath system, retrieval loop system, Coons dilator, needle, wire guide and three-way stopcock.



www.cookmedical.com

COOK MEDICAL INCORPORATED
P.O. Box 4195, Bloomington, IN 47402-4195 U.S.A.
Phone: 812.339.2235, Toll Free: 800.457.4500, Toll Free Fax: 800.554.8335

COOK (CANADA) INC.
111 Sandiford Drive, Stouffville, Ontario, L4A 7X5 CANADA
Phone: 905.640.7110, Toll Free: 800.668.0300

WILLIAM A. COOK AUSTRALIA PTY. LTD.
95 Brandl Street, Brisbane Technology Park, Eight Mile Plains
Brisbane, QLD 4113 AUSTRALIA, Phone: +61 7 3841 1188

WILLIAM COOK EUROPE ApS
Sandet 6, DK-4632, Bjaeverskov, DENMARK, Phone: +45 56 86 86 86

© COOK 2010 PIV-BM-CVCFNDS-EN-201002