



SILICONE BALLOON HSG CATHETER

Achieve maximum visualization of
the uterus and fallopian tubes.



CONTROL

Small diameter
allows easier passage
through the cervix.

VISIBILITY

Size of balloon enables
enhanced view during
evaluation.

VERSATILITY

Catheter can be used
to perform both SHG
and HSG procedures.



SILICONE BALLOON HSG CATHETER



The Cook® Silicone Balloon HSG Catheter is used to access the uterine cavity for saline infusion sonohysterography or to inject contrast medium into the uterine cavity and fallopian tubes for evaluation of tubal patency. Sold in boxes of 10.

- Silicone construction; contains no latex.
- Available with or without sideport.
- Includes an appropriate syringe.

Cook Medical—Developing simple solutions for women’s health.

Order Number	Reference Number	Fr	Length cm	Balloon Volume mL	Syringe Volume mL
Closed Tip with Sideport					
G17472	J-CHSG-503001	5.0	30	1.0	1.0
G17471	J-CHSG-703001	7.0	30	1.5	3.0
Open Tip					
G17147	J-CHSG-503000	5.0	30	1.0	1.0
G17870	J-CHSG-703000	7.0	30	1.5	3.0

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

Customer Service Centers

- AORTIC INTERVENTION
- CRITICAL CARE
- ENDOSCOPY
- INTERVENTIONAL RADIOLOGY
- LEAD MANAGEMENT
- PERIPHERAL INTERVENTION
- SURGERY
- UROLOGY
- WOMEN'S HEALTH**

EMEA: EDI - www.cookmedical.com/edi.do
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Poland: +48 223060159, pl.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335,
 E-mail: orders@cookmedical.com
Australia: ABN: 51 141259118
 Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283,
 E-mail: wacustserv@cookust.com.au



www.cookmedical.com