

Access and orientation.



Cotton CannulaTome II®

PRECURVED DOUBLE LUMEN SPHINCTEROTOME

With its excellent orientation and smooth DomeTip®, the Cannulatome II pc allows you to access the ductal system more easily and with potentially reduced trauma. Plus, on our Protector™ model, we've added insulation on the proximal portion of the cutting wire to reduce the incidence of contact between the cutting wire and endoscope, as well as an on-off sidearm adaptor to conveniently convert to triple-lumen functionality.



www.cookmedical.com



Cotton CannulaTome II*

PRECURVED DOUBLE LUMEN SPHINCTEROTOME

Used for endoscopic cannulation of the ductal system and for sphincterotomy. The device is supplied sterile and is disposable - intended for single use only.

- Precurved DomeTip provides easier access with potentially less trauma*
- Patented 3-D forming wire preserves optimal tip shape during shipping and storage, assuring superb orientation during the procedure
- Protector™ insulation on proximal portion of cutting wire decreases incidence of cutting wire contact with scope, reducing the risk of damage to the cutting wire and scope (available on CT-25M-P)
- On-off sidearm adaptor conveniently delivers triple-lumen functionality to double-lumen device (available on CT-25M-P)

Order Number	Reference Part Number	Tip	Cutting Wire Length mm	Cutting Wire Type	Wire Guide Diameter* in	Color Connection
Cotton CannulaTome® II PC Protector™						
G49196	CT-25M-P	DomeTip	25	Monofilament	.035	Purple
Cotton CannulaTome® II PC						
G22241	CT-25	DomeTip	25	Braided	.035	Purple
G22537	CT-25M	DomeTip	25	Monofilament	.035	Purple
CannulaTome II® PC						
G22240	CT-20	DomeTip	20	Braided	.035	Purple
G22536	CT-20M	DomeTip	20	Monofilament	.035	Purple
G25801	CT-30	DomeTip	30	Braided	.035	Purple
G22538	CT-30M	DomeTip	30	Monofilament	.035	Purple

Minimum accessory channel 2.8 mm. 6 Fr catheter *Recommend Tracer Metro® Wire, sold separately. Active cord available separately.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

*Data on file at Cook Medical.

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: orders@cookmedical.com

Australia:
 Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283
 E-mail: cau.custserv@cookmedical.com

