

Seeing isn't just believing,  
it's everything.



## EchoTip<sup>®</sup> Ultra

HD ENDOSCOPIC ULTRASOUND NEEDLE




[www.cookmedical.com](http://www.cookmedical.com)

# Available in high definition.


In today's complex, exacting EUS procedures, needle visualization is more important than ever. The better the needle visualization, the more confident you can be.

That's why Cook created the EchoTip® Ultra with HDFNA (High Definition Fine Needle Aspiration) - the most advanced EUS needle yet.

In fact, tests show that this new, unique echogenic HD pattern - makes it be up to 3 times brighter than other needles currently on the market.\*



Highly echogenic needle allows you to more adequately assure the needle tip is within the target, which may improve your cytology yields



• Nitinol stylet for stability and flexibility


• A naturally contoured handle delivers precise needle control and stability

• "Zero" reference mark ensures complete needle retraction within the sheath

**Up to 3 times brighter than other needles currently on the market.\***

# High definition in the palm of your hands.

Complementing the superb echogenicity at its distal end, the handle of the EchoTip Ultra with HDFNA has the advanced features of the original EchoTip Ultra that allow you to confidently perform any EUS procedure referred to you.






• Safety lock ring easily slides and locks at the desired extension

• Integrated sheath adjuster offers compatibility with a full range of echoendoscopes

• Metal Luer Lock for secure anchoring to the echoendoscope

# Expanding your view

No matter what needle size your ultrasound procedure requires, a quality image is key to successful outcomes. That's why Cook now offers the EchoTip Ultra with HDFNA in 19 and 25 gage to go along with its 22 gage needle. Three different platforms - all in high definition



# Product Specifications

Used to sample targeted submucosal gastrointestinal lesions through the accessory channel of an ultrasound endoscope.

Supplied sterile and is disposable - intended for single use only.

Order Number	Reference Part Number	Needle Gage	Sheath Size Fr	Needle Extension cm	Stylet Tip	Dimple Pattern
G31518	ECHO-1-22	22	5.2	0-8	2 mm Ball Tip	High Definition
G31521	ECHO-3-22	22	5.2	0-8	2 mm Ball Tip	High Definition
G31520	ECHO-19	19	5.2-4.2	0-8	2 mm Ball Tip	High Definition
G31519	ECHO-25	25	5.2	0-8	2 mm Ball Tip	High Definition

*Minimum accessory channel 2.0 mm. Adjustable sheath length. 5 cc and 10 cc vacuum locking syringe included.*

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

## Customer Service

**EMEA:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)  
Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)  
Belgium: +32 27001633, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)  
Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)  
Finland: +358 972519996, [fi.orders@cookmedical.com](mailto:fi.orders@cookmedical.com)  
France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)  
Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)  
Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)  
Norway: +47 23162968, [no.orders@cookmedical.com](mailto:no.orders@cookmedical.com)  
Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)  
Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)  
Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

[www.cookmedical.com](http://www.cookmedical.com)

**Americas:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335  
E-mail: [orders@cookmedical.com](mailto:orders@cookmedical.com)

**Australia:**  
Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283  
E-mail: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)

