



FALLOPIAN TUBE CATHETERIZATION SET

Catheterize, inject and evaluate fallopian tubes with one device.



COMPREHENSIVE

All components for tubal catheterization are included in one set.

SPECIFIC

Designed for a distinct part of the anatomy and tailor-made for its indication.

EFFECTIVE

Microcatheter enhances tubal access through the guide catheter.



FALLOPIAN TUBE CATHETERIZATION SET



The Fallopian Tube Catheterization Set is used for selective catheterization of the proximal fallopian tube(s), injection of contrast medium, and evaluation of tubal patency. The coaxial catheter assembly advances into position through a device such as the Intrauterine Access Balloon or other uterine access device.

- Set includes two wire guides for tubal catheterization.
- Designed to be placed through any Cook uterine access device with a 6 Fr or larger working channel.
- Blue microcatheter is intended to enhance tubal access as it moves inside the gray torque control catheter.

Cook Medical—Developing simple solutions for women’s health.

Order Number	Reference Number	Torque Control Catheter Fr	Injection Catheter Fr	Cope Mandril Wire Guide inch / mm
G10991	FTC-550	5.5	3.0	.015/.38
G10992	FTC-550-NT	5.5	3.0	.018/.46

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

Customer Service Centers

- AORTIC INTERVENTION
- CRITICAL CARE
- ENDOSCOPY
- INTERVENTIONAL RADIOLOGY
- LEAD MANAGEMENT
- PERIPHERAL INTERVENTION
- SURGERY
- UROLOGY
- WOMEN'S HEALTH**

EMEA: EDI - www.cookmedical.com/edi.do
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Poland: +48 223060159, pl.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335,
 E-mail: orders@cookmedical.com
Australia: ABN: 51 141259118
 Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283,
 E-mail: wcacustserv@cookaust.com.au



www.cookmedical.com