


**COOK**  
MEDICAL

# Fusion<sup>®</sup> Zilver<sup>®</sup>

**BILIARY SELF-EXPANDING STENT**




Palliation is all about maintaining patency.



# Fusion® Zilver®

## BILIARY SELF-EXPANDING STENT



Used in palliation of malignant neoplasms of the biliary tree.

**WARNING:** The safety and effectiveness of this device for use in the vascular system have not been established.

Supplied sterile and is disposable - intended for single use only.

- Self-expanding, non-shortening nitinol construction provides maximum flexibility without kinking or compressing
- Zilver "Z" configuration applies uniform radial strength evenly to the biliary duct
- Multiple gold markers on each end for precise placement
- Compatible with all wire guide lengths
- Preloaded 7.0 FR Reliant Introducer offers an integrated hub for injection with Tuohy-Borst adapter to prevent inadvertent deployment

Order Number	Reference Part Number	Stent Diameter mm	Stent Length cm	Introducer French size/ Length Fr/cm	Wire Guide Diameter inch	Minimum Accessory Channel mm
G31555	FS-ZILBS-6-4	6	4	7/208	.035	4.2
G31556	FS-ZILBS-8-4	8	4	7/208	.035	4.2
G31557	FS-ZILBS-10-4	10	4	7/208	.035	4.2
G31558	FS-ZILBS-6-6	6	6	7/208	.035	4.2
G31559	FS-ZILBS-8-6	8	6	7/208	.035	4.2
G31560	FS-ZILBS-10-6	10	6	7/208	.035	4.2
G31561	FS-ZILBS-6-8	6	8	7/208	.035	4.2
G31562	FS-ZILBS-8-8	8	8	7/208	.035	4.2
G31563	FS-ZILBS-10-8	10	8	7/208	.035	4.2

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

### Customer Service Centers

**EMEA:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
 Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)  
 Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)  
 Belgium: +32 27001633, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)  
 Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)  
 France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)  
 Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)  
 Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)  
 Norway: +47 23162968, [no.orders@cookmedical.com](mailto:no.orders@cookmedical.com)  
 Poland: +48 223060159, [pl.orders@cookmedical.com](mailto:pl.orders@cookmedical.com)  
 Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)  
 Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)  
 Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

[www.cookmedical.com](http://www.cookmedical.com)

**Americas:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335,  
 E-mail: [orders@cookmedical.com](mailto:orders@cookmedical.com)

**Australia:**  
 Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283,  
 E-mail: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)

- AORTIC INTERVENTION
- CRITICAL CARE
- ENDOSCOPY
- INTERVENTIONAL RADIOLOGY
- LEAD MANAGEMENT
- PERIPHERAL INTERVENTION
- SURGERY
- UROLOGY
- WOMEN'S HEALTH

