


# Dilate with strength and accuracy.



## Hercules<sup>®</sup>

3 STAGE WIRE GUIDED BALLOON

To effectively manage strictures in the esophagus, pylorus, duodenum and colon, you need a strong balloon that can be accurately staged. When you choose the Cook Hercules, you get both. Thanks to the exclusive P.E.T. Flex™ technology, this balloon balances exceptional strength and staging accuracy. And now, for clinician preference, the Hercules includes a preloaded wire guide. The wire-guided Hercules is yet another addition to Cook's full line of innovative endoscopic accessories designed to meet all your patients' needs.



[www.cookmedical.com](http://www.cookmedical.com)



## Hercules®

### 3 STAGE WIRE GUIDED BALLOON

Used to endoscopically dilate strictures of the gastrointestinal tract including the esophagus, pylorus, duodenum and colon.

Supplied sterile and is disposable - intended for single use only.

- Rapid deflation completes procedure quickly\*
- Two radiopaque markers identify body of balloon for excellent visibility of balloon positioning
- Simplified packaging enhances procedural efficiency
- Glow-in-the-dark size and inflation information on packaging and catheter tags for easy identification
- 2.8 mm minimum accessory channel for all scope manufacturers

Order Number	Reference Part Number	Inflated Balloon Diameter mm	Inflated Balloon Diameter Fr	Balloon Length cm	ATM Pressure	PSI Pressure
G51799	HBD-W-8-9-10	8-9-10	24-27-30	5.5	4-6-8	60-90-120
G51795	HBD-W-10-11-12	10-11-12	30-33-36	5.5	2-4-6	30-60-90
G51796	HBD-W-12-13.5-15	12-13.5-15	36-40.5-45	5.5	2-4-6	30-60-90
G51797	HBD-W-15-16.5-18	15-16.5-18	45-49.5-54	5.5	2-4-6	30-60-90
G51798	HBD-W-18-19-20	18-19-20	54-57-60	5.5	2-4-6	30-60-90

**Note:** Catheter size is 8 Fr. Catheter length 240 cm. Minimum accessory channel 2.8 mm.

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

\*Based on using a 60cc syringe. Data on file at Cook Medical.

### Customer Service

**EMEA:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
 Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)  
 Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)  
 Belgium: +32 27001633, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)  
 Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)  
 Finland: +358 972519996, [fi.orders@cookmedical.com](mailto:fi.orders@cookmedical.com)  
 France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)  
 Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)  
 Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)  
 Norway: +47 23162968, [no.orders@cookmedical.com](mailto:no.orders@cookmedical.com)  
 Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)  
 Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)  
 Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

[www.cookmedical.com](http://www.cookmedical.com)

**Americas:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335  
 E-mail: [orders@cookmedical.com](mailto:orders@cookmedical.com)

**Australia:**  
 Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283  
 E-mail: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)

