

# Use the MReye<sup>®</sup> coil as an anchor in large and high-flow vessels.



## MReye<sup>®</sup> EMBOLIZATION COIL

### High radial force

The high radial force of the MReye provides strong purchase and stability in large or high-flow vessels.

### Extra-long fibers

Fibers on the highly visible Inconel<sup>®</sup> MReye coil allow you to achieve rapid occlusion in high-flow vessels.

Inconel is a registered trademark of the Huntington Alloys Corporation.

### Multiple configurations

Designed for .035 inch and .038 inch systems, the MReye coil comes in a variety of lengths (2-15 cm) and diameters (2-25 mm) so that you can occlude large or high-flow vessels.



[www.cookmedical.com](http://www.cookmedical.com)

## MREYE EMBOLIZATION COIL


Order Number	Reference Part Number	Recommended Catheter ID/ End Hole Diameter inch	Extended Embolus Length cm	Coiled Embolus Diameter mm	Approximate Number of Loops
<b>.035 inch</b>					
G42406	IMWCE-35-2-3	.035	2	3	2.1
G47900	IMWCE-35-3-2	.035	3	2	4.8
G42408	IMWCE-35-3-3	.035	3	3	3.1
G42409	IMWCE-35-3-4	.035	3	4	2.4
G42410	IMWCE-35-3-5	.035	3	5	1.9
G42414	IMWCE-35-4-3	.035	4	3	4.2
G42419	IMWCE-35-5-3	.035	5	3	5.3
G36405	IMWCE-35-5-5	.035	5	5	3.1
G20291	IMWCE-35-5-6	.035	5	6	2.7
G19907	IMWCE-35-5-7	.035	5	7	2.3
G42423	IMWCE-35-5-8	.035	5	8	2.0
G19908	IMWCE-35-5-10	.035	5	10	1.6
G19957	IMWCE-35-8-5	.035	8	5	5.1
G20948	IMWCE-35-8-8	.035	8	8	3.2
G20207	IMWCE-35-8-10	.035	8	10	2.5
G20452	IMWCE-35-10-8	.035	10	8	4.0
G20035	IMWCE-35-10-10	.035	10	10	3.2
G47887	IMWCE-35-12-10	.035	12	10	3.8
G47894	IMWCE-35-14-8	.035	14	8	5.6
G47890	IMWCE-35-14-12	.035	14	12	3.7
G19929	IMWCE-35-15-10	.035	15	10	4.8
G47895	IMWCE-35-15-12	.035	15	12	4.0
G20055	IMWCE-35-15-15	.035	15	15	3.2
G19949	IMWCE-35-15-20	.035	15	20	2.4

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

## MREYE EMBOLIZATION COIL

Order Number	Reference Part Number	Recommended Catheter ID/ End Hole Diameter inch	Extended Embolus Length cm	Coiled Embolus Diameter mm	Approximate Number of Loops
<b>.038 inch</b>					
G42448	IMWCE-38-2-3	.038	2	3	2.1
G42451	IMWCE-38-3-5	.038	3	5	1.9
G42454	IMWCE-38-4-3	.038	4	3	4.2
G42458	IMWCE-38-5-3	.038	5	3	5.3
G42460	IMWCE-38-5-5	.038	5	5	3.2
G42463	IMWCE-38-5-8	.038	5	8	2.0
G42465	IMWCE-38-5-10	.038	5	10	1.6
G42466	IMWCE-38-5-12	.038	5	12	1.3
G42467	IMWCE-38-5-15	.038	5	15	1.1
G42472	IMWCE-38-8-5	.038	8	5	5.1
G42475	IMWCE-38-8-10	.038	8	10	2.5
G42477	IMWCE-38-8-15	.038	8	15	1.7
G19920	IMWCE-38-10-10	.038	10	10	3.2
G47914	IMWCE-38-12-10	.038	12	10	3.8
G47917	IMWCE-38-12-14	.038	12	14	2.7
G19886	IMWCE-38-15-10	.038	15	10	4.8
G42486	IMWCE-38-15-12	.038	15	12	4.0
G42487	IMWCE-38-15-15	.038	15	15	3.2
G19887	IMWCE-38-15-20	.038	15	20	2.4
G42489	IMWCE-38-15-25	.038	15	25	1.9

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.



## MReye®

EMBOLIZATION COIL

- MReye coils are used for arterial and venous embolization in the peripheral vasculature.
- MReye coils maintain higher radial force than platinum varieties.
- The radiodensity of the coils is higher than stainless steel. This makes the coils more visible under fluoroscopy.
- MReye coils are MR conditional up to 3T, SAR of 2.0 W/kg.
- The coils come in .035 and .038 inch diameters.
- MReye coils are available in extended embolus lengths of 2-15 cm and coiled embolus lengths of 3-25 mm.
- The approximate number of loops on the coils ranges from 1.1 to 5.6.

### Customer Service

**EMEA:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)

Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)

Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)

Belgium: +32 27001633, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)

Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)

Finland: +358 972519996, [fi.orders@cookmedical.com](mailto:fi.orders@cookmedical.com)

France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)

Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)

Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)

Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)

Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)

Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)

Norway: +47 23162968, [no.orders@cookmedical.com](mailto:no.orders@cookmedical.com)

Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)

Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)

Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)

Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)

Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)

United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

[www.cookmedical.com](http://www.cookmedical.com)

**Americas:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: [orders@cookmedical.com](mailto:orders@cookmedical.com)

**Australia:**

Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283

E-mail: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)

