Take complete control of your PEG placement procedures.



Percutaneous Endoscopic Gastrostomy System

From the moment you open the Cook PEG kit, you'll notice the difference. The components are conveniently arranged in the sequential order in which they will be used, making the procedure flow as smoothly and efficiently as the nutrients you are delivering to your patients. The Cook PEG kit - available in push or pull methods - delivers the efficiency you need and the nutrients your patients need.



Percutaneous Endoscopic Gastrostomy System

Used for percutaneous endoscopic placement to provide enteral nutrition to patients requiring nutritional support.

Supplied sterile and is disposable - intended for single use only.

- Kit components, including hemostat and tube clamp, are conveniently packaged for sequential use during the procedure for increased efficiency
- Clear feeding tube made of flexible, high-grade medical silicone is optimal length to give physician control throughout placement
- Sharper scissors provide for a better cut during the placement procedure
- Reference "X mark" to specify where to cut

Order Number	Reference Part Number	Feeding Tube Fr	Feeding Tube Inner Diameter mm	Insertion Wire Diameter inch	Insertion Wire Length cm
G22637	FLOW-20-PUSH-S	20	4.47	.035	260
G22638	FLOW-20-PULL-S	20	4.47	n/a	150
G22635	PEG-24-PUSH-S	24	5.5	.035	260
G22636	PEG-24-PULL-S	24	5.5	n/a	150
G31543	FLOW-20-PUSH-SSK-S*	20	4.47	.035	260
G56403	FLOW-20-PULL-SSK-S-T*	20	4.47	n/a	150
G31548	PEG-24-PUSH-SSK-S*	24	5.5	.035	260
G56402	PEG-24-PULL-SSK-S-T*	24	5.5	n/a	150
G22762	FLOW-20-PUSH-I-S**	20	4.47	.035	260
G22761	FLOW-20-PULL-I-S**	20	4.47	n/a	150
G22758	PEG-24-PUSH-I-S**	24	5.5	.035	260
G22757	PEG-24-PULL-I-S**	24	5.5	n/a	150

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

KIT CONTENTS:	1 Silicone Feeding Tube	1 Universal Adapter	1 Patient Care Manual	1 25 gage x 5/8" Needle**
1 Fenestrated Drape	5 Gauze pads, 4 x 4"	1 Bolus Adapter	1 #11 Scalpel*	1 ChloraPrep® Swabs (USA Only)**
1 150 cm Looped Insertion	1 Scissors	1 Twist Lock	1 Trocar Needle*	1 Povidone Ointment (USA Only)**
Wire (Pull Only)	1 Hemostat	1 Cable Tie	2 Water Soluble Lubricant	
1 260 cm Wire guide (Push Only)	1 Bolster	Bolus Feeding Adapter – Male (International Only)	(USA Only)**	(USA Only)*,**
1 Cold Snare	1 Tube Clamp		1 22 gage x 1½" Needle**	1 Syringe*,**

^{*}Safety components available within Safety Sharp Kits (-SSK) **Kits without pharmaceuticals (-I Kits)

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do Distributors: +353 61239240, ssc.distributors@cookmedical.com Austria: +43 179567121, oe.orders@cookmedical.com Belgium: +32 27001633, be.orders@cookmedical.com Denmark: +45 38487607, da.orders@cookmedical.com Finland: +358 972519996, fi.orders@cookmedical.com France: +33 171230269, fr.orders@cookmedical.com Germany: +49 6950072804, de.orders@cookmedical.com Hungary: +36 17779199, hu.orders@cookmedical.com Ireland: +353 61239252, ie.orders@cookmedical.com Italy: +39 0269682853, it.orders@cookmedical.com Netherlands: +31 202013367, nl.orders@cookmedical.com Norway: +47 23162968, no.orders@cookmedical.com Spain: +34 912702691, es.orders@cookmedical.com Sweden: +46 858769468, se.orders@cookmedical.com Switzerland - French: +41 448009609, fr.orders@cookmedical.com Switzerland - Italian: +41 448009609, it.orders@cookmedical.com Switzerland - German: +41 448009609, de.orders@cookmedical.com United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: orders@cookmedical.com

Australia:

Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283

E-mail: cau.custserv@cookmedical.com