

Curved Introducer

## **CONVENIENT REMOVAL**

Dual knobs allow the sheath to be easily peeled away after insertion is completed.

## TRUSTED DESIGN

Cook Medical invented peel-away sheath technology that has now become the industry standard.

## COMPLETE SETS

Both curved and straight sets are available. In addition to introducer sheath, each set includes radiopaque dilator (with special van Andel taper on curved sets); .038 inch, double-ended stainless steel wire guide; 18 gage access needle; and 12 mL syringe.

Pass pacemaker leads or catheters <u>smoothly into the ven</u>ous system.





The curved Peel-Away® Sheath Introducer Set allows pacemaker leads or catheters to be delivered to a desired location distal to the superior vena cava/innominate juncture. The 44 cm set is especially useful for coronary sinus leads, and the other curved sets are especially useful for reinsertions.

Supplied sterile in peel-open packages. Intended for one-time use.

| Global<br>Product<br>Number | Order<br>Number          | Sheath<br>ID<br>Fr | Length<br>cm |
|-----------------------------|--------------------------|--------------------|--------------|
| Curved Introd               | ucer Set                 |                    |              |
| G31000                      | PLVN-7.0-38-24-VAD5-CVI  | 7.0                | 24           |
| G29841                      | PLVN-9.0-38-24-VAD5-CVI  | 9.0                | 24           |
| G29838                      | PLVN-10.0-38-24-VAD5-CVI | 10.0               | 24           |
| G29839                      | PLVN-12.0-38-24-VAD5-CVI | 12.0               | 24           |
| G29840                      | PLVN-14.0-38-24-VAD5-CVI | 14.0               | 24           |
| G29690                      | PLVN-9.0-38-44-VAD5-CVI  | 9.0                | 44           |
| Straight Intro              | ducer Set                |                    |              |
| G04456                      | C-PLI-7.0-38             | 7.0                | 15.5         |
| G01922                      | C-PLI-8.0-38             | 8.0                | 15.5         |
| G04502                      | C-PLI-9.0-38             | 9.0                | 15.5         |
| G04495                      | C-PLI-10.0-38            | 10.0               | 15.5         |
| G04496                      | C-PLI-11.0-38            | 11.0               | 15.5         |



COOK MEDICAL INCORPORATED PO. Box 4195, Bloomington, IN 47402-4195 U.S.A. Phone: 812.339.2235, Toll Free: 800.457.4500, Toll Free Fax: 800.554.8335 COOK (CANADA) INC. 111 Sandiford Drive, Stouffville, Ontario, L4A 7X5 CANADA Phone: 905.640.7110, Toll Free: 800.668.0300

www.cookmedical.com

WILLIAM A. COOK AUSTRALIA PTY. LTD. 95 Brandl Street, Brisbane Technology Park, Eight Mile Plains Brisbane, QLD 4113 AUSTRALIA, Phone: +61 7 3841 1188

WILLIAM COOK EUROPE ApS Sandet 6, DK-4632, Bjaeverskov, DENMARK, Phone: +45 56 86 86 86

LM-BOUSM-PADS-EN-201006

© COOK 2010

| AORTIC       |  |
|--------------|--|
| INTERVENTION |  |

CRITICAL ENDOSCOPY

CARE

INTERVENTIONAL RADIOLOGY

LEAD MANAGEMENT

PERIPHERAL INTERVENTION

SURGERY UROLOGY

WOMEN'S HEALTH