



# Royal Flush® Plus

BEACON® TIP HIGH-FLOW CATHETERS



Enhance your view during catheter manipulation and injection.

#### HIGH FLOW RATES

High-strength nylon construction with larger inner diameters facilitate optimal flow rates.

#### SUPERIOR VISUALIZATION

Cook's Beacon® enhanced radiopaque tip technology ensures excellent visibility during catheter manipulation.

#### ATRAUMATIC TIP DESIGN

Soft tip configurations enhance transition to the wire guide.



[www.cookmedical.com](http://www.cookmedical.com)

**COOK MEDICAL INCORPORATED**

P.O. Box 4195, Bloomington, IN 47402-4195 U.S.A.  
Phone: 812.339.2235, Toll Free: 800.457.4500, Toll Free Fax: 800.554.8335

**COOK (CANADA) INC.**

111 Sandiford Drive, Stouffville, Ontario, L4A 7X5 CANADA  
Phone: 905.640.7110, Toll Free: 800.668.0300

**WILLIAM A. COOK AUSTRALIA PTY. LTD.**

95 Brandl Street, Brisbane Technology Park, Eight Mile Plains  
Brisbane, QLD 4113 AUSTRALIA, Phone: +61 7 3841 1188

**WILLIAM COOK EUROPE ApS**

Sandet 6, DK-4632, Bjaeverskov, DENMARK, Phone: +45 56 86 86 86

© COOK 2010      PIL-2M-HNRCDS-EN-201007

AORTIC  
INTERVENTION

CRITICAL  
CARE

ENDOSCOPY

INTERVENTIONAL  
RADIOLOGY

LEAD  
MANAGEMENT

PERIPHERAL  
INTERVENTION

SURGERY

UROLOGY


WOMEN'S  
HEALTH



# Royal Flush® Plus

BEACON® TIP HIGH-FLOW CATHETERS

Global Product Number	Order Number	Fr	Accepts Wire Guide Diameter inch	Length cm	Sideports	Maximum Flow Rate mL/sec	Maximum Pressure kg/cm <sup>2</sup>	Maximum Pressure psi
<b>Verrill Contralateral Flush</b>								
G12394	HNR4.0-35-65-P-8S-VCF	4.0	.035	65	8	19	84	1,200
G12999	HNR4.0-35-90-P-8S-VCF	4.0	.035	90	8	17	84	1,200
G12430	HNR4.0-35-100-P-6S-VCF	4.0	.035	100	6	15	84	1,200
G12129	HNR5.0-35-65-P-8S-VCF	5.0	.035	65	8	32	84	1,200
G12799	HNR5.0-35-90-P-8S-VCF	5.0	.035	90	8	29	84	1,200
G31403	HNR5.0-35-100-P-8S-VCF	5.0	.035	100	8	27	84	1,200
<b>Pigtail Flush</b>								
G09406	HNR4.0-35-70-P-10S-PIG	4.0	.035	70	10	20	84	1,200
G09407	HNR4.0-35-90-P-10S-PIG	4.0	.035	90	10	17	84	1,200
G11278	HNR4.0-35-100-P-10S-PIG	4.0	.035	100	10	15	84	1,200
G11423	HNR4.0-35-110-P-10S-PIG	4.0	.035	110	10	14	84	1,200
G08832	HNR5.0-35-70-P-10S-PIG	5.0	.035	70	10	32	84	1,200
G09408	HNR5.0-35-90-P-10S-PIG	5.0	.035	90	10	29	84	1,200
G09059	HNR5.0-35-100-P-10S-PIG	5.0	.035	100	10	27	84	1,200
G11265	HNR5.0-35-110-P-10S-PIG	5.0	.035	110	10	26	84	1,200
G11272	HNR5.0-38-70-P-10S-PIG	5.0	.038	70	10	32	84	1,200
G11275	HNR5.0-38-90-P-10S-PIG	5.0	.038	90	10	29	84	1,200
G09285	HNR5.0-38-100-P-10S-PIG	5.0	.038	100	10	27	84	1,200
G11269	HNR5.0-38-110-P-10S-PIG	5.0	.038	110	10	26	84	1,200
<b>Pediatric Pigtail Flush</b>								
G11263	HNR4.0-35-70-P-4S-PED	4.0	.035	70	4	17	84	1,200
G11728	HNR4.0-35-90-P-6S-PED	4.0	.035	90	6	17	84	1,200
G11281	HNR5.0-35-70-P-4S-PED	5.0	.035	70	4	29	84	1,200
G11729	HNR5.0-35-90-P-6S-PED	5.0	.035	90	6	29	84	1,200






# Royal Flush® Plus

BEACON® TIP HIGH-FLOW CATHETERS



Global Product Number	Order Number	Fr	Accepts Wire Guide Diameter inch	Length cm	Sideports	Maximum Flow Rate mL/sec	Maximum Pressure kg/cm <sup>2</sup>	Maximum Pressure psi
<b>Straight Flush</b>								
G09409	HNR4.0-35-70-P-10S-0	4.0	.035	70	10	20	84	1,200
G09410	HNR4.0-35-90-P-10S-0	4.0	.035	90	10	17	84	1,200
G09111	HNR4.0-35-100-P-10S-0	4.0	.035	100	10	15	84	1,200
G09411	HNR5.0-35-70-P-10S-0	5.0	.035	70	10	32	84	1,200
G09412	HNR5.0-35-90-P-10S-0	5.0	.035	90	10	29	84	1,200
G09697	HNR5.0-35-100-P-10S-0	5.0	.035	100	10	27	84	1,200
G11241	HNR5.0-35-110-P-10S-0	5.0	.035	110	10	26	84	1,200
G11270	HNR5.0-38-70-P-10S-0	5.0	.038	70	10	32	84	1,200
G11273	HNR5.0-38-90-P-10S-0	5.0	.038	90	10	29	84	1,200
G11267	HNR5.0-38-100-P-10S-0	5.0	.038	100	10	27	84	1,200
G09461	HNR5.0-38-110-P-10S-0	5.0	.038	110	10	26	84	1,200



<b>Centra-Flo™ Pigtail Flush</b>								
G11279	HNR4.0-35-70-P-10S-CFP	4.0	.035	70	10	20	84	1,200
G11151	HNR4.0-35-90-P-10S-CFP	4.0	.035	90	10	17	84	1,200
G11422	HNR4.0-35-100-P-10S-CFP	4.0	.035	100	10	15	84	1,200
G11262	HNR5.0-35-70-P-10S-CFP	5.0	.035	70	10	32	84	1,200
G11266	HNR5.0-35-90-P-10S-CFP	5.0	.035	90	10	29	84	1,200
G11264	HNR5.0-35-100-P-10S-CFP	5.0	.035	100	10	27	84	1,200
G11271	HNR5.0-38-70-P-10S-CFP	5.0	.038	70	10	32	84	1,200
G11274	HNR5.0-38-90-P-10S-CFP	5.0	.038	90	10	29	84	1,200
G11268	HNR5.0-38-100-P-10S-CFP	5.0	.038	100	10	27	84	1,200