




SMALL GAGE OVUM ASPIRATION NEEDLES



RETRIEVE



Enhance accuracy in oocyte collection.




SMALL GAGE OVUM ASPIRATION NEEDLES

Collecting immature oocytes requires precision and control.

The Small Gage Ova-Stiff™ Ovum Aspiration Needle from Cook Medical offers a finely tuned bevel tip specifically designed for effectiveness in small antral follicles. Highly visible under ultrasonic imaging, this tip enables accuracy, and the stiff cannula gives the clinician further control of placement. As a result, more viable oocytes can be collected.

Cook Medical—Specializing in positive patient outcomes.



CookARTlab.com

Features and benefits of ovum aspiration needles:



Small Gage Ova-Stiff

- Extra-stiff cannula for precise placement of tip in follicle.
- Standard RWH noncoring bevel prevents tissue from entering lumen.
- Translucent FEP vacuum line. Flushing is achieved by attaching the flushing cannula to a syringe and inserting the cannula into the flared end of the aspiration line at the silicone stopper.
- EchoTip® enhances visualization of needle tip when used with ultrasonic imaging equipment.



Immature Ovum Aspiration Needle

- EchoTip extends around the bevel, improving visibility and placement of the needle.
- 17 gage, 30 cm guide needle used to puncture vaginal wall, adding rigidity to the 20 gage, 35 cm aspiration needle.
- New stopper design accepts the tip of a Luer slip syringe, facilitating follicle flushing.



Small Gage Ova-Stiff™

OVUM ASPIRATION NEEDLE

Global Product Number	Order Number	gage	Length cm	Aspiration Line cm	Bevel	Handle Color
G20943	K-OPS-6035-RWH-B-ET	18	35	60	B	pink
G19715	K-OPS-6035-RWH-ET	18	35	60	RWH	pink
G26182	K-OPS-7030-RWH-ET	19	30	60	RWH	brown
G26638	K-OPS-7035-RWH-ET	19	35	60	RWH	brown
G26631	K-OPS-8035-RWH-ET	20	35	60	RWH	yellow
G26056	K-OPS-8035-RWH-B-ET	20	35	60	B	yellow
G50735	K-OPS-2130*	21	30	75	RWH	green
G50736	K-OPS-2135*	21	35	75	RWH	green

*Fitted with 50 cm vacuum line



IMMATURE OVUM ASPIRATION NEEDLE

Global Product Number	Order Number	Guide Needle gage	Guide Needle Length cm	Aspiration Needle gage	Aspiration Needle Length cm	Aspiration Tubing Length cm
G38609	K-IOPS-2035-1730	17	30	20	35	60

Fertility specialists around the world can rely on Cook Medical's expertly crafted assisted reproduction products. Each one is engineered through collaborative research to meet the evolving and exacting needs of this delicate process and optimize outcomes for patients and their families.



RETRIEVE
CULTURE
TRANSFER



www.cookmedical.com

COOK MEDICAL INCORPORATED

P.O. Box 4195, Bloomington, IN 47402-4195 USA

Phone: 812.339.2235, Toll Free: 800.457.4500, Toll Free Fax: 800.554.8335

COOK IRELAND LTD.

O'Halloran Road, National Technology Park,

Limerick, IRELAND, Phone: +353 61 334440, Fax: +353 61 334441

COOK (CANADA) INC.

111 Sandiford Drive, Stouffville, Ontario, L4A 7X5 CANADA

Phone: 905.640.7110, Toll Free: 800.668.0300, Fax: 905.640.6804

WILLIAM A. COOK AUSTRALIA PTY. LTD.

95 Brandl Street, Brisbane Technology Park, Eight Mile Plains

Brisbane, QLD 4113 AUSTRALIA, Phone: +61 7 3841 1188

© COOK 2009

WH-BM-SGNB-EN-200906

AORTIC
INTERVENTION

CARDIOLOGY

CRITICAL
CARE

ENDOSCOPY

INTERVENTIONAL
RADIOLOGY

PERIPHERAL
INTERVENTION

SURGERY

UROLOGY

WOMEN'S
HEALTH