


# Confidently navigate difficult anatomy.




## Tracer Metro<sup>®</sup>

WIRE GUIDE

Effective treatment is your ultimate destination and the Tracer Metro can take you there. The seamless hydrophilic coated radiopaque tip navigates smoothly through difficult anatomy. Vivid reference markings help you monitor wire movement and facilitate stent selection. Uniquely designed for biliary access, the Tracer Metro gives you an efficient balance of flexibility and visibility.



[www.cookmedical.com](http://www.cookmedical.com)



## Tracer Metro®

### WIRE GUIDE

Used to assist in cannulation and in bridging difficult strictures during ERCP.

Supplied sterile and is disposable - intended for single use only.

Placement of the Tracer Metro allows for continual access to the desired duct and facilitates exchange of compatible wire guided products.

- 5 cm seamless hydrophilic coated radiopaque tip provides clear and smooth navigation with enhanced radiopacity
- Kink-resistant core provides enhanced maneuverability
- Numerical indicators at 1 cm increments allow you to reference stricture length beginning distally 6 to 24 cm

Order Number	Reference Part Number	Wire Guide Length cm	Wire Guide Diameter inch	Radiopaque Tip	Color Connection
G25263	MET-25-480	480	.025	Straight	Blue
G25264	MET-25-480-A	480	.025	Angled	Blue
G25212	MET-35-480	480	.035	Straight	Purple
G25213	MET-35-480-A	480	.035	Angled	Purple

5 cm Aqua Coat Tip, radiopaque tip. May be left in place during sphincterotomy. **Caution:** Not for use with metal tip devices.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

### Customer Service

**EMEA:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
 Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)  
 Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)  
 Belgium: +32 27001633, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)  
 Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)  
 Finland: +358 972519996, [fi.orders@cookmedical.com](mailto:fi.orders@cookmedical.com)  
 France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)  
 Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)  
 Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)  
 Norway: +47 23162968, [no.orders@cookmedical.com](mailto:no.orders@cookmedical.com)  
 Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)  
 Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)  
 Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

[www.cookmedical.com](http://www.cookmedical.com)

**Americas:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335  
 E-mail: [orders@cookmedical.com](mailto:orders@cookmedical.com)

**Australia:**  
 Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283  
 E-mail: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)

