


# Aspirate with enhanced visualization.



## ECHOTIP® DISPOSABLE AMNIOCENTESIS NEEDLE

### Visibility

EchoTip® technology allows the needle to be seen more clearly under ultrasound.

### Stability


Strong, sturdy needles provide rigidity for better handling.

### Convenience

Multiple gages and lengths accommodate different body types.



[www.cookmedical.com](http://www.cookmedical.com)



## ECHOTIP® DISPOSABLE AMNIOCENTESIS NEEDLE

The EchoTip® Disposable Amniocentesis Needle is used to aspirate fluid from the amniotic sac. Designed exclusively for use with ultrasound, this needle comes in a variety of gages and lengths.

- Needle design provides enhanced visualization of the needle tip with ultrasonic imaging equipment.
- Available with or without a sideport.
- Distal sideport aids in fluid aspiration.

Cook Medical—Developing simple solutions for women’s health.

Order Number	Reference Part Number	gage	Length cm
Without Sideport			
G16283	J-DAN-200901	20	9
G16287	J-DAN-201201	20	12
G16289	J-DAN-201501	20	15
G16291	J-DAN-202001	20	20
G16297	J-DAN-211201	21	12
G16299	J-DAN-211501	21	15
G35572	J-DAN-212001	21	20
G16303	J-DAN-220901	22	9
G16307	J-DAN-221201	22	12
G16309	J-DAN-221501	22	15
G35573	J-DAN-222001	22	20

Order Number	Reference Part Number	gage	Length cm
With Sideport			
G16284	J-DAN-200902	20	9
G16288	J-DAN-201202	20	12
G16290	J-DAN-201502	20	15
G16294	J-DAN-210902	21	9
G16298	J-DAN-211202	21	12
G16300	J-DAN-211502	21	15
G16304	J-DAN-220902	22	9
G16308	J-DAN-221202	22	12
G16310	J-DAN-221502	22	15

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

### Customer Service

**EMEA:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
 Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)  
 Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)  
 Belgium: +32 27001633, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)  
 Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)  
 Finland: +358 972519996, [fi.orders@cookmedical.com](mailto:fi.orders@cookmedical.com)  
 France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)  
 Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)  
 Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)  
 Norway: +47 23162968, [no.orders@cookmedical.com](mailto:no.orders@cookmedical.com)  
 Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)  
 Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)  
 Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

[www.cookmedical.com](http://www.cookmedical.com)

**Americas:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335  
 E-mail: [orders@cookmedical.com](mailto:orders@cookmedical.com)

**Australia:**  
 Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283  
 E-mail: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)

