



Lunderquist[®]

EXTRA-STIFF WIRE GUIDE



Guide with optimal support.

An Accessory for Endovascular Procedures

Straight Curved Double Curved Extended Curved



The Lunderquist® Extra-Stiff Wire Guide is designed to provide maximum support for navigation of tortuous anatomy during endovascular procedures. With four highly flexible wire-guide tip configurations available, Lunderquist provides compatibility and unrivaled support.

Support

A PTFE-coated, stainless steel mandril is the key to Lunderquist providing maximum support and low surface friction during complex diagnostic and interventional procedures, including endovascular aortic repair.

Safety

The smooth transition from the stiff mandril to the flexible gold tip aids atraumatic vessel navigation. A radiopaque gold inner coil allows precise tip visibility and location during wire-guide manipulation.

Compatibility

Lunderquist is compatible with a wide variety of anatomies. Various lengths and four tip configurations are available, including double-curved and extended-curve tip designs that mimic the anatomy.

Order Number	Reference Part Number	Diameter inch	Length cm	Flexible Tip Length cm	Tip Curve Radius mm	Coil Length cm
Straight Exchange - PTFE-Coated Stainless Steel						
G31453	TSMG-35-260-LES	.035	260	4	straight	15
G45352	TSMG-35-300-LES	.035	300	4	straight	15
Curved Exchange - PTFE-Coated Stainless Steel						
G45353	TSCMG-35-260-7-LES	.035	260	4	7.5	15
G45407	TSCMG-35-300-7-LES	.035	300	4	7.5	15
Order Number	Reference Part Number	Diameter inch	Length cm	Flexible Tip Length cm	Primary/Secondary Curve Radius mm	Coil Length cm
Double Curved Exchange - PTFE-Coated Stainless Steel						
G45208	TSCMG-35-260-LESDC	.035	260	4	75/15	15
G45408	TSCMG-35-300-LESDC	.035	300	4	75/15	15
Extended Curve Exchange - PTFE-Coated Stainless Steel						
G53967	TSCMG-35-260-E-LESDC	.035	260	4	55/15	15
G53971	TSCMG-35-300-E-LESDC	.035	300	4	55/15	15

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

Customer Service Centers

AORTIC INTERVENTION
CRITICAL CARE
ENDOSCOPY
INTERVENTIONAL RADIOLOGY
LEAD MANAGEMENT
PERIPHERAL INTERVENTION
SURGERY
UROLOGY
WOMEN'S HEALTH

EMEA: EDI - www.cookmedical.com/edi.do
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001633, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 972519996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
E-mail: orders@cookmedical.com

Australia:
Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283
E-mail: cau.custserv@cookmedical.com