





Zenith® Fenestrated

AAA ENDOVASCULAR GRAFT



Transcend: Precision

Experience patient-specific precision.





Zenith® Fenestrated

AAA ENDOVASCULAR GRAFT

It's about precision.


The Zenith® Fenestrated AAA Endovascular Graft is **the first graft specifically designed for precise positioning** in patients with infrarenal necks as short as 4 mm. Now with the friction-reducing and kink-resistant features of the Flexor® Introducer Sheath, Zenith Fenestrated makes treatment of more AAA patients even easier.

Going beyond. That's what it means to Transcend. That's the essence of Zenith.


precision Design

During the manufacturing of Zenith Fenestrated, precisely placed fenestrations and scallops are positioned to accommodate an individual patient's unique renal artery and SMA anatomy.


This patient-specific composition allows uncompromised positioning of the graft's proximal sealing stent to achieve an optimal seal in patients with short necks.



precision Placement



- ① **Flexor Introducer Sheath**
Kink-resistant technology with low-friction hydrophilic coating.
- ② **AQ® Hydrophilic Coating**
Dramatically reduces friction to enhance deliverability.
- ③ **Diameter-Reducing Ties**
Constrain graft diameter during deployment to facilitate placement and visceral vessel cannulation.
- ④ **Fenestrations and Scallops**
Equipped with gold markers for added visibility.
- ⑤ **Trigger Wires**
Facilitate staged deployment to ensure precise placement.



To learn more about the unrivaled migration resistance provided by ARC Technology™, please visit www.cookmedical.com/arc.

AORTIC INTERVENTION
CRITICAL CARE
ENDOSCOPY
INTERVENTIONAL RADIOLOGY
LEAD MANAGEMENT
PERIPHERAL INTERVENTION
SURGERY
UROLOGY
WOMEN'S HEALTH



Customer Service Centers

EMEA: EDI - www.cookmedical.com/edi.do
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001633, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Poland: +48 223060159, pl.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
E-mail: orders@cookmedical.com

Australia:
Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283
E-mail: cau.custserv@cookmedical.com

