



# CENTESIS & DRAINAGE



Discover a comprehensive solution  
to meet your drainage needs.



# Contents

## **Pneumothorax Catheters**

- 3 Pneumothorax Sets
- 4 Wayne Pneumothorax Sets
- 5 Cook® Emergency Pneumothorax Set
- 6 Simple Pneumothorax Aspiration Accessory Set

## **Multipurpose Catheters**

- 7 Thal-Quick Chest Tube Sets - Single Lumen
- 8 Thal-Quick Chest Tube Trays - Single Lumen
- 9 Thal-Quick Chest Tube Set and Tray - Double Lumen
- 10 Thal-Quick Chest Tube Adapter
- 11 Fuhrman Pleural/Pneumopericardial Drainage Set

## **Pericardiocentesis Catheters**

- 12 Pericardiocentesis Set
- 13 Lock Pericardiocentesis Set
- 14 Pigtail Pericardiocentesis Tray
- 15 Teitel Pediatric Pericardiocentesis Set

## **Peritoneal Lavage Catheters**

- 16 Peritoneal Lavage Set and Tray
- 17 Educational Materials
- 18 Resources and Acknowledgments

## Pneumothorax Sets

Used to treat spontaneous and tension pneumothorax.



### FEATURES & BENEFITS

- Trocar placement.
- Adapter on tubing can be connected to a suction apparatus or to the Cook Chest Drain Valve.
- Radiopaque catheter enhances X-ray visualization.

### COMPONENTS

- Needle
- Radiopaque catheter
- Cook Chest Drain Valve
- One-way stopcock
- Connecting tube

Intended for one-time use.

Global Product Number	Order Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Catheter Material	Needle gage	Needle Length cm
<b>Pediatric</b>							
G03465	C-TPTP-100	6.3	19	6	PTFE	21	12
<b>Standard</b>							
G03301	C-TPT-100	9.0	29	20	PTFE	18	20
G03302	C-TPT-200	10.2	29	20	polyurethane	18	20

## Wayne Pneumothorax Sets

Used for the relief of simple, spontaneous, iatrogenic and tension pneumothorax.



### FEATURES & BENEFITS

- Seldinger or trocar placement.
- Percutaneous (Seldinger) placement facilitates controlled, minimally invasive catheter introduction.
- Enlarged sideports enhance drainage of air.
- Radiopaque catheter enhances X-ray visualization.

### COMPONENTS

#### Seldinger Set

- EchoTip® echogenic introducer needle
- 9 cm introducer needle
- Wire guide
- Dilator
- Pigtail catheter
- Cook Chest Drain Valve

#### Trocar Set

- Needle obturator
- Pigtail catheter
- Cook Chest Drain Valve

Intended for one-time use.

Global Product Number	Order Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Needle gage	Needle Length cm
<b>Seldinger Set*</b>						
G11182	C-UTPT-1400-WAYNE-112497	14.0	29	19	18	5, 9
<b>Trocar Sets*</b>						
G09922	C-UTPT-1020-WAYNE	10.2	29	20	19	16.5
G07910	C-UTPT-1400-WAYNE	14.0	29	20	19	16.5

\*Sets also include scalpel, three-way stopcock and connecting tube.

## Cook® Emergency Pneumothorax Set

Used for emergency relief and temporary management of suspected tension pneumothorax.



### FEATURES & BENEFITS

- Catheter-over-needle placement.
- Reinforced catheter prevents collapse after placement.
- Cook Chest Drain Valve permits continued re-expansion of the lung by natural forces, preventing inadvertent backflow of air into the pleural space.

### COMPONENTS

- Catheter introducer needle
- Cook Chest Drain Valve
- One-way stopcock
- Connecting tube
- Catheter fixation device
- Syringe
- Transparent tape
- Alcohol prep
- Chlorhexidine

Intended for one-time use.

Global Product Number	Order Number	Catheter Fr	Catheter Length cm	Needle gage
G13722	C-TPTS-8.5-6.0-FSNS	8.5	6	15

## Simple Pneumothorax Aspiration Accessory Set

Used in conjunction with pneumothorax aspiration catheters if incomplete expansion of the lung persists after mechanical aspiration of air.



### FEATURES & BENEFITS

- Cook Chest Drain Valve permits continued re-expansion of the lung by natural forces, preventing inadvertent backflow of air into the pleural space.

### COMPONENTS

- Cook Chest Drain Valve
- Connecting tube

Intended for one-time use.

Global Product Number	Order Number	Connecting Tube Fr	Connecting Tube Length cm
G05297	C-CASP-A-FORD	14.0	30

## Thal-Quick Chest Tube Sets - Single Lumen

Used for pleural fluid drainage.



### FEATURES & BENEFITS

- Percutaneous (Seldinger) placement facilitates controlled, minimally invasive catheter introduction.
- EchoTip® echogenic introducer needle enhances needle tip visibility under ultrasound imaging.
- Radiopaque catheter enhances X-ray visualization.
- Marked dilator facilitates controlled insertion of chest tube into the pleural cavity.

### SET COMPONENTS

- Scalpel
- EchoTip echogenic introducer needle
- Amplatz Extra-Stiff Wire Guide
- Marked dilator
- Chest tube inserter
- Thal-Quick Chest Tube

Intended for one-time use.

Global Product Number	Order Number	Fr	Length cm	Sideports
<b>Sets</b>				
G06786	C-TQTS-800	8.0	18	3
G06787	C-TQTS-1000	10.0	20	3
G05161	C-TQTS-1200	12.0	22	3
G05777	C-TQTS-1400	14.0	41	4
G04891	C-TQTS-1600	16.0	41	4
G05453	C-TQTS-1800	18.0	41	4
G05454	C-TQTS-2000	20.0	41	4
G04137	C-TQTS-2400	24.0	41	4
G04139	C-TQTS-2800	28.0	41	4
G07243	C-TQTS-3200	32.0	41	4
G07241	C-TQTS-3600	36.0	41	4



## Thal-Quick Chest Tube Trays - Single Lumen

Used for pleural fluid drainage.



### TRAY COMPONENTS

- Scalpel
- EchoTip® echogenic introducer needle
- Amplatz Extra-Stiff Wire Guide
- Dilators
- Chest tube inserter
- Thal-Quick Chest Tube
- 25 gage needle
- 22 gage needle
- Lidocaine
- Chlorhexidine
- Fenestrated drape
- Gauze sponges
- Disposable syringes
- Suture with needle
- Povidone-iodine ointment
- CSR wrap
- Prep tray
- Needle holder cup

Intended for one-time use.

Global Product Number	Order Number	Fr	Length cm	Sideports
<b>Trays</b>				
G06791	C-TQTSY-800	8.0	18	3
G06885	C-TQTSY-1000	10.0	20	3
G05464	C-TQTSY-1200	12.0	22	3
G07090	C-TQTSY-1400	14.0	41	4
G05462	C-TQTSY-1600	16.0	41	4
G05461	C-TQTSY-1800	18.0	41	4
G05463	C-TQTSY-2000	20.0	41	4
G04220	C-TQTSY-2400	24.0	41	4
G04226	C-TQTSY-2800	28.0	41	4
G07242	C-TQTSY-3200	32.0	41	4
G07244	C-TQTSY-3600	36.0	41	4

## Thal-Quick Chest Tube Set and Tray - Double Lumen

Used for pleural fluid drainage.



### FEATURES & BENEFITS

- Additional lumen available for aspiration or infusion.
- Sideports are positioned at the distal end of the chest tube.
- Radiopaque catheter enhances X-ray visualization.

### COMPONENTS

- Scalpel
- Introducer needle
- Amplatz Extra-Stiff Wire Guide
- Marked dilators
- Chest tube inserter
- Thal-Quick Chest Tube
- Needle holder cup
- Multipurpose tubing adapter

Intended for one-time use.

Global Product Number	Order Number	Fr	Length cm	Sideports
<b>Set</b>				
G36224	C-DLCTS-2400	24.0	40	4 large, 4 small
<b>Tray*</b>				
G36225	C-DLCTSY-2400	24.0	40	4 large, 4 small

\* Tray also includes the following items: 25 gage needle, 22 gage needle, lidocaine, chlorhexidine, fenestrated drape, gauze sponges, disposable syringes, suture with needle, povidone-iodine ointment, CSR wrap, prep tray and needle holder cup.

## Thal-Quick Chest Tube Adapter

Used to provide separate sampling or infusion port within Thal-Quick Chest Tubes.



### FEATURES & BENEFITS

- Accepts 24.0-36.0 Fr devices

### COMPONENTS

- Thal-Quick Chest Tube Adapter
- Three-way stopcock

Intended for one-time use.

Global Product Number	Order Number	Accepts Fr
G07467	C-TQF-28.0	24.0-36.0

## Fuhrman Pleural/Pneumopericardial Drainage Set

Used for evacuation of air from the pericardial sac or to drain air or fluid from the pleural space.



### FEATURES & BENEFITS

- Percutaneous (Seldinger) placement facilitates controlled, minimally invasive catheter introduction.
- Sideports positioned within catheter pigtail.
- Centimeter markings on pigtail catheter.
- Radiopaque catheter enhances X-ray visualization.

### COMPONENTS

- Introducer needle
- Wire guide
- Dilator
- Radiopaque pigtail catheter
- Three-way stopcock
- Multipurpose tubing adapter

Intended for one-time use.

Global Product Number	Order Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Wire Guide Diameter inch	Wire Guide Length cm	Wire Guide Tip Configuration	Needle gage	Needle Length cm
G03974	C-PPD-850	8.5	15	6	.028	40	2 mm J	18	7

## Pericardiocentesis Set

Used to drain fluid from the pericardial sac.



### FEATURES & BENEFITS

- Percutaneous (Seldinger) placement facilitates controlled, minimally invasive catheter introduction.
- Stainless steel hub of introducer sheath needle permits constant electrocardiographic monitoring via alligator clip cable.
- Radiopaque catheter enhances X-ray visualization.

### COMPONENTS

- Sheath needle
- Wire guide
- Stiffening cannula
- Radiopaque catheter
- Alligator clip cable
- Connecting tube, 7.5 cm long
- Connecting tube, 5 cm long

Intended for one-time use.

Global Product Number	Order Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Stiffening Cannula gage	Stiffening Cannula Length cm	Wire Guide Diameter inch	Wire Guide Length cm	Wire Guide Tip Configuration	Needle gage	Needle Length cm
G03282	C-PCS-850	8.5	22	4	18	24	.038	70	3 mm J	18	15

## Lock Pericardiocentesis Set

Used to drain fluid from the pericardial sac.



### FEATURES & BENEFITS

- Percutaneous (Seldinger) placement facilitates controlled, minimally invasive catheter introduction.
- Radiopaque catheter enhances X-ray visualization.

### COMPONENTS

- Introducer needles
- Wire guide
- Dilator
- Radiopaque pigtail catheter
- Multipurpose tubing adapter

Intended for one-time use.

Global Product Number	Order Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Wire Guide Diameter inch	Wire Guide Length cm	Wire Guide Tip Configuration	Needle gage	Needle Length cm
G02433	C-PCS-830-LOCK	8.3	40	6	.038	70	3 mm J	18	7, 12

## Pigtail Pericardiocentesis Tray

Used to drain fluid from the pericardial sac.



### FEATURES & BENEFITS

- Alligator clip cable is attached to the hub of the needle for constant electrocardiographic monitoring during needle insertion.
- Includes tubing and specimen tubes.

### COMPONENTS

- Introducer needle
- Wire guide
- Dilators
- Radiopaque pigtail catheter
- Three-way stopcock
- Alligator clip cable
- Multipurpose tubing adapter
- 15 mL conical specimen tubes
- Connecting tube

Intended for one-time use.

Global Product Number	Order Number*	Catheter Fr	Catheter Length cm	Catheter Sideports	Wire Guide Diameter inch	Wire Guide Length cm	Wire Guide Tip Configuration	Needle gage	Needle Length cm
G05252	C-PCSY-830-AT	8.3	40	6	.038	70	straight and 3 mm J	18	8

\* Tray also includes the following items: 25 gage needle, 22 gage needle, lidocaine, chlorhexidine, fenestrated drape, gauze sponges, disposable syringes, disposable scalpel, suture with needle, povidone-iodine ointment, CSR wrap, prep tray and needle holder cup.

## Teitel Pediatric Pericardiocentesis Set

Used to drain fluid from the pericardial sac.



### FEATURES & BENEFITS

- Percutaneous (Seldinger) placement facilitates controlled, minimally invasive catheter introduction.
- Stainless steel hub of introducer needle permits constant electrocardiographic monitoring via alligator clip cable.
- Radiopaque catheter enhances X-ray visualization.

### COMPONENTS

- Introducer needle
- Wire guide
- Dilator
- Radiopaque pigtail catheter
- Alligator clip cable
- Connecting tube, 7.5 cm long
- Connecting tube, 5 cm long

Intended for one-time use.

Global Product Number	Order Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Wire Guide Diameter inch	Wire Guide Length cm	Wire Guide Tip Configuration	Needle gage	Needle Length cm
G05251	C-PCS-500-TTL	5.0	30	6	.028	60	3 mm J	19	6



## Peritoneal Lavage Set and Tray

Used to infuse and withdraw rinsing solutions in the peritoneal cavity.



### FEATURES & BENEFITS

- Percutaneous (Seldinger) placement facilitates controlled, minimally invasive catheter introduction.
- 90 sideports on catheter enhance flow of lavage fluid.
- Radiopaque catheter enhances X-ray visualization.

### COMPONENTS

- Introducer needle
- Straight safety wire guide
- Radiopaque catheter

Intended for one-time use.


Global Product Number	Order Number	Catheter Fr	Catheter Length cm	Catheter Sideports	Wire Guide Diameter inch	Wire Guide Length cm	Needle gage	Needle Length cm
<b>Set</b>								
G05267	C-PLS-901	9.0	20	90	.038	40	18	7
<b>Tray*</b>								
G05268	C-PLSY-901	9.0	20	90	.038	40	18	7

\*Tray also includes the following items: 25 gage needle, 22 gage needle, lidocaine, chlorhexidine, fenestrated drape, gauze sponges, disposable syringes, disposable scalpel, suture with needle, povidone-iodine ointment, CSR wrap, prep tray and needle holder cup.

## Centesis & Drainage Products

CD-ROM

CD-BD-CAD-ML-200904



## Resources

Sargent EN, Turner AF. Emergency treatment of pneumothorax. A simple catheter technique for use in the radiology department. *Am J Roentgenol Radium Ther Nucl Med.* 1970;109(3):531-535.

Lock JE, Bass JL, Kulik TJ, et al. Chronic percutaneous pericardial drainage with modified pigtail catheters in children. *Am J Cardiol.* 1984;53(8):1179-1182.

Fuhrman BP, Landrum BG, Ferrara TB, et al. Pleural drainage using modified pigtail catheters. *Crit Care Med.* 1986;14(6):575-576.

Tilkian AG, Daily EK. Pericardiocentesis and drainage. In: Tilkian AG, Daily EK, eds. *Cardiovascular Procedures: Diagnostic Techniques and Therapeutic Procedures.* St. Louis, MO: C.V. Mosby Company, 1986.

Thal AP, Quick KL. A guided chest tube for safe thoracostomy. *Surg Gynecol Obstet.* 1988;167(6):517.

Vallee P, Sullivan M, Richardson H, et al. Sequential treatment of a simple pneumothorax. *Ann Emerg Med.* 1988;17(9):936-942.

Conces DJ Jr, Tarver RD, Gray WC, et al. Treatment of pneumothoraces utilizing small caliber chest tubes. *Chest.* 1988;94(1):55-57.

Wood B, Dubik M. A new device for pleural drainage in newborn infants. *Pediatrics.* 1995;96(5, pt 1):955-956.

Ahmed MY, Silver P, Nimkoff L, et al. The needle-wire-dilator technique for the insertion of chest tubes in pediatric patients. *Pediatr Emerg Care.* 1995;11(4):252-254.

Higginson I. Equipment reviews: emergency pneumothorax set. *Pre-Hosp Immed Care.* 1998;2:104-105.

Gammie JS, Banks MC, Fuhrman CR, et al. The pigtail catheter for pleural drainage: a less invasive alternative to tube thoracostomy. *JSL.* 1999;3(1):57-61.

## Acknowledgments

CF Landers, MD. Division of Pulmonary and Critical Care Medicine, Sharp Memorial Hospital, San Diego, California.

DF Teitel, MD. Department of Pediatric Cardiology, University of California at San Francisco, San Francisco, California.



[www.cookmedical.com](http://www.cookmedical.com)

**COOK MEDICAL INCORPORATED**

P.O. Box 4195, Bloomington, IN 47402-4195 U.S.A.  
Phone: 812.339.2235, Toll Free: 800.457.4500, Toll Free Fax: 800.554.8335

**COOK (CANADA) INC.**

111 Sandiford Drive, Stouffville, Ontario, L4A 7X5 CANADA  
Phone: 905.640.7110, Toll Free: 800.668.0300

**WILLIAM A. COOK AUSTRALIA PTY. LTD.**

95 Brandl Street, Brisbane Technology Park, Eight Mile Plains  
Brisbane, QLD 4113 AUSTRALIA, Phone: +61 7 3841 1188

**WILLIAM COOK EUROPE ApS**

Sandet 6, DK-4632, Bjaeverskov, DENMARK, Phone: +45 56 86 86 86

© COOK 2010

CC-BUSM-CADB-EN-201003

AORTIC  
INTERVENTION

CRITICAL  
CARE

ENDOSCOPY

INTERVENTIONAL  
RADIOLOGY

LEAD  
MANAGEMENT

PERIPHERAL  
INTERVENTION

SURGERY

UROLOGY

WOMEN'S  
HEALTH