


Choose the only SFA stent
with a proven drug effect.



Zilver® PTX®
DRUG-ELUTING PERIPHERAL STENT




www.cookmedical.com

Proven drug effect.

Established Technology

Zilver® PTX® is based on the next-generation Zilver Flex stent platform, which is known for its balance of flexibility, fracture resistance, and radial force. The combination of drug-elution and a strong yet flexible scaffold is the **standard of care for SFA lesions.**



**Next-
Generation
Durable Design**


Fracture Rate of
0.9%

Superior Results

A randomised trial demonstrated the **clear drug effect** of Zilver PTX. At three years, Zilver PTX achieved a 79.6% patency rate, compared with 56.3% for provisional bare-metal stents.¹

3 Year Drug Effect - Patency (PSVR < 2.0)

(Zilver PTX vs. Provisional BMS)



Group	3-Year Restenosis ¹	
	Rate	Reduction
Provisional Zilver PTX	20.4%	53%
Provisional Bare Zilver	43.7%	

1. Dake M. The Zilver PTX randomized trial of paclitaxel-eluting stents for femoropopliteal disease: 3-year results. Presented at: Vascular Interventional Advances (VIVA) 2012; October 9-12, 2012; Las Vegas, Nevada.

Longer-Lasting Treatment

Zilver PTX provides sustained effectiveness. With Zilver PTX, you can achieve longer-lasting outcomes that are superior to the outcomes achieved by bare-metal stents.

CASE 1

54-year-old male diabetic


Preprocedure



Postprocedure



4.9 Years later



CASE 2

63-year-old male diabetic


Preprocedure



Postprocedure



4.8 Years later



Images courtesy of Dr. Fabrizio Fanelli,
University of Rome, Italy.

20,000 +
patients,

6 continents,
and the

1 and only
drug-eluting
SFA stent.

Zilver® PTX®

DRUG-ELUTING PERIPHERAL STENT


Order Number	Reference Part Number	Wire Guide Diameter inch	Stent Diameter mm	Stent Length mm	Minimum Sheath Fr
80 cm Introducer Length					
G55284	ZIV6-35-80-5.0-40-PTX	.035	5	40	6
G55285	ZIV6-35-80-5.0-60-PTX	.035	5	60	6
G55286	ZIV6-35-80-5.0-80-PTX	.035	5	80	6
G55287	ZIV6-35-80-5.0-100-PTX	.035	5	100	6
G55291	ZIV6-35-80-6.0-40-PTX	.035	6	40	6
G55292	ZIV6-35-80-6.0-60-PTX	.035	6	60	6
G55293	ZIV6-35-80-6.0-80-PTX	.035	6	80	6
G55294	ZIV6-35-80-6.0-100-PTX	.035	6	100	6
G55301	ZIV6-35-80-7.0-40-PTX	.035	7	40	6
G55302	ZIV6-35-80-7.0-60-PTX	.035	7	60	6
G55303	ZIV6-35-80-7.0-80-PTX	.035	7	80	6
G55304	ZIV6-35-80-7.0-100-PTX	.035	7	100	6
G55298	ZIV6-35-80-8.0-40-PTX	.035	8	40	6
G55272	ZIV6-35-80-8.0-60-PTX	.035	8	60	6
G55273	ZIV6-35-80-8.0-80-PTX	.035	8	80	6
G55274	ZIV6-35-80-8.0-100-PTX	.035	8	100	6

Order Number	Reference Part Number	Wire Guide Diameter inch	Stent Diameter mm	Stent Length mm	Minimum Sheath Fr
125 cm Introducer Length					
G51333	ZIV6-35-125-5.0-40-PTX	.035	5	40	6
G51334	ZIV6-35-125-5.0-60-PTX	.035	5	60	6
G51335	ZIV6-35-125-5.0-80-PTX	.035	5	80	6
G51336	ZIV6-35-125-5.0-100-PTX	.035	5	100	6
G51341	ZIV6-35-125-6.0-40-PTX	.035	6	40	6
G51342	ZIV6-35-125-6.0-60-PTX	.035	6	60	6
G51343	ZIV6-35-125-6.0-80-PTX	.035	6	80	6
G51344	ZIV6-35-125-6.0-100-PTX	.035	6	100	6
G51349	ZIV6-35-125-7.0-40-PTX	.035	7	40	6
G51350	ZIV6-35-125-7.0-60-PTX	.035	7	60	6
G51351	ZIV6-35-125-7.0-80-PTX	.035	7	80	6
G51352	ZIV6-35-125-7.0-100-PTX	.035	7	100	6
G51356	ZIV6-35-125-8.0-40-PTX	.035	8	40	6
G51357	ZIV6-35-125-8.0-60-PTX	.035	8	60	6
G51358	ZIV6-35-125-8.0-80-PTX	.035	8	80	6
G51359	ZIV6-35-125-8.0-100-PTX	.035	8	100	6

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Order Number Key

ZIV6-35-80-6.0-100-PTX



Customer Service

EMEA: EDI - www.cookmedical.com/edi.do
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: orders@cookmedical.com

Australia:
 Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283
 E-mail: cau.custserv@cookmedical.com

