


# Tip flexibility: Changing with the times

One of your procedural goals will never change: simple, reliable access during ERCP. One thing that has changed is tip flexibility. Over time, wire guides have become more and more flexible. At Cook, we didn't stop until we made a wire guide tip that is measurably one of the most flexible on the market: the Acrobat Calibrated Tip wire guide.



Timeline of Innovation

# Product Specifications

Used to assist in cannulation of the biliary and pancreatic ducts and to aid in bridging difficult strictures during ERCP.

Supplied sterile and is disposable – intended for single use only.

Order Number	Reference Part Number	Wire Diameter inch	Wire Length cm	Radiopaque Tip cm	Hydrophilic Coating cm
G34264	ACRO-35-205	.035	205	4	25
G34265	ACRO-35-260	.035	260	4	25
G34266	ACRO-35-450	.035	450	4	25
G32758	ACRO-25-260	.025	260	4	25
G32759	ACRO-25-450	.025	450	4	25

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

## Acrobat Preloaded Sphincterotomes

Used for cannulation of the ductal system and for sphincterotomy. If preloaded, also aids in bridging difficult strictures during ERCP.

Supplied sterile and is disposable – intended for single use only.

Order Number	Reference Part Number	Catheter Size Fr	Tip Shape	Wire Length cm	Wire Diameter inch	Minimum Accessory Channel mm
G34601	FS-OMNI-ACRO-35-205	7	DomeTip	205	.035	3.2
G34693	FS-OMNI-ACRO-35-260	7	DomeTip	260	.035	3.2
G34711	FS-OMNI-ACRO-35-450	7	DomeTip	450	.035	3.2
G34780	DASH-ACRO-35-260	6.5	DomeTip	260	.035	2.8
G34778	DASH-ACRO-35-450	6.5	DomeTip	450	.035	2.8
G34779	DASH-ACRO-25-450	6	DomeTip	450	.025	2.8

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

### Customer Service

**EMEA:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
 Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)  
 Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)  
 Belgium: +32 27001633, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)  
 Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)  
 Finland: +358 972519996, [fi.orders@cookmedical.com](mailto:fi.orders@cookmedical.com)  
 France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)  
 Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)  
 Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)  
 Norway: +47 23162968, [no.orders@cookmedical.com](mailto:no.orders@cookmedical.com)  
 Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)  
 Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)  
 Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

[www.cookmedical.com](http://www.cookmedical.com)

**Americas:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335  
 E-mail: [orders@cookmedical.com](mailto:orders@cookmedical.com)

**Australia:**  
 Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283  
 E-mail: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)



Performance you can measure.



Acrobat™

CALIBRATED TIP WIRE GUIDE



[www.cookmedical.com](http://www.cookmedical.com)




# Interventional Radiology access now for ERCP

The same performance that sets the standard for wire guide access in interventional radiology is now available for ERCP. We tailored interventional radiology technology specifically for ERCP wire guide access. With one of the most flexible tips on the market, Acrobat gives you the same Cook-pioneered technology that interventionalists rely on everyday for difficult selective and superselective vascular access.



# Performance you can measure

The only sure way to prove wire guide tip flexibility is to measure it. That's what we did, and Acrobat performed superbly. It is so flexible that it required the least amount of gram-force to initiate an alpha loop 1 cm from the distal tip. So now, whether you're cannulating or traversing difficult strictures, the Acrobat Calibrated Tip wire guide can be your everyday wire guide that performs like a specialty wire.



\*Data on file at Cook Medical.



Back-loading indicator allows smooth transition to next device of choice



Spring-coil design for one of the most flexible tips on the market, helping you gain access and navigate difficult strictures



Endoscopic markings indicate wire guide movement and allows measurement using the Mark V System



Fluoroscopic images courtesy of Dr. Stan Branch  
Duke University Hospital, Durham, NC.

Endoscopic image courtesy of Dr. Adam Slivka  
UPMC Presbyterian Hospital, Pittsburgh, PA.