


Redefine access.



Aprima™ Access


NONVASCULAR INTRODUCER SET



www.cookmedical.com

Redefine access with a set that is engineered to make every aspect of your procedure go smoothly.

Transitionless-Tip™



The Aprima Access set's Transitionless-Tip design requires less insertion force than standard access sets¹ and virtually eliminates hang-ups during entry, which helps provide seamless access and reduce the risk of patient trauma.

Full-Length Radiopacity




The entire shaft and distal tip—not just one small band—are radiopaque to maximize fluoroscopic visibility during placement.

Ease of Placement



The set includes an EchoTip® echogenic access needle for optimal ultrasound visibility, a wire guide, and a hydrophilic-coated coaxial introducer sheath to ease every step of the placement process.


¹Data on file.



Aprima™ Access

NONVASCULAR INTRODUCER SET

- Used for single-puncture percutaneous access to facilitate placement of a working wire guide for drainage procedures.
- The Transitionless-Tip design, hydrophilic coating, and set components work together to ease placement.
- The complete sets include a coaxial introducer and a choice of two types of EchoTip echogenic needles and two types of wire guides in order to suit almost any physician preference.



Order Number	Reference Part Number	Sheath ID/OD Fr	Needle gage/Length cm	Needle Tip Configuration	Wire Guide Material
G34253	NSSW-4.0-18-NPAS-100-HC-SST*	4.0/6.0	-	-	-
G52962	NPAS-100-HC-U-SST	4.0/6.0	22/15	Chiba	stainless steel
G52963	NPAS-100-HC-NT-U-SST	4.0/6.0	22/15	Chiba	nitinol
G52964	NPAS-104-HC-U-SST	4.0/6.0	21/15	trocac	stainless steel
G52965	NPAS-104-HC-NT-U-SST	4.0/6.0	21/15	trocac	nitinol
G52966	NPAS-104-HC-NT-JWG-U-SST**	4.0/6.0	21/15	trocac	nitinol
G52967	NPAS-105-HC-NT-U-SST	4.0/6.0	22/20	Chiba	nitinol
G25516	NPAS-108-HC-NT-U-SST	4.0/6.0	21/15	Chiba	nitinol
G25517	NPAS-120-HC-NT-U-SST	4.0/6.0	20/15	Chiba	nitinol

*This set does not include a needle or a wire guide.

**This set includes a .038 inch, 100 cm heavy-duty wire guide with a curved tip.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.



Customer Service

EMEA: EDI - www.cookmedical.com/edi.do
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001633, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 972519996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
E-mail: orders@cookmedical.com

Australia:
Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283
E-mail: cau.custserv@cookmedical.com

