

Bridge difficult strictures with efficiency.



D.A.S.H.[®] ERCP CATHETER

When facing challenging ERCP procedures, you need a catheter that performs efficiently. That's why the D.A.S.H. ERCP Catheter incorporates an integrated hub and a Tuohy-Borst adapter for simultaneous wire guide access and contrast injection. You can also choose a version that includes a .025" Tracer Metro[®] Direct[™] Wire Guide preloaded in a racetrack conveniently attached to the hub.



www.cookmedical.com



D.A.S.H.[®]

ERCP CATHETER

This device is used for endoscopic cannulation of the ductal system.

The device is supplied sterile and is disposable - intended for single use only.

- Tuohy-Borst adapter provides wire guide access while simultaneously allowing contrast injection
- Offered preloaded with Tracer Metro[®] Direct[™] wire guide, non-preloaded accepts .025" wire guide
- Wire guide racetrack connection to catheter hub (in DASH-T-480 only) simplifies procedure for assistant, reducing procedure time
- 5.5 FR tapered to 4.5 FR catheter compatible with 2.8 mm endoscope
- Mark V System provides endoscopic visualization at 5, 10 and 15 cm from the distal tip of the catheter signifying cannulation depth
- Radiopaque tip enhances fluoroscopic visualization

Order Number	Reference Part Number	French Size Fr	Catheter Tip Size Fr	Catheter Length cm	Wire Guide Diameter inch	Color Connection
G25514	DASH-T	5.5	4.5	200	.025*	Blue
G25515	DASH-T-480	5.5	4.5	200	Preloaded .025	Blue

*Recommend Tracer Metro[®] Direct[™] Wire Guide, sold separately. Minimum accessory channel 2.8 mm.

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

Customer Service Centers

EMEA: EDI - www.cookmedical.com/edi.do
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: orders@cookmedical.com

Australia:
 Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283
 E-mail: cau.custserv@cookmedical.com



www.cookmedical.com