

# Take your staging procedures to the next level.



## EchoTip® Ultra

HD ENDOBRONCHIAL ULTRASOUND NEEDLE

In today's complex, exacting EBUS procedures, needle visualization is more important than ever. That's why Cook, the leader in endoscopic ultrasound, created the EchoTip® Ultra Endobronchial High Definition Ultrasound Needle—a uniquely designed needle that can potentially enhance your diagnostic and staging capabilities. The needle's high echogenicity allows you to more adequately assure the needle tip is within the target, which may improve cytological yields. The ergonomic handle and scope adapter can also optimize procedural performance, stability and security.



[www.cookmedical.com](http://www.cookmedical.com)



Olympus



Pentax

## EchoTip® Ultra

HD ENDOBRONCHIAL ULTRASOUND NEEDLE

**Olympus** - Used to sample targeted submucosal and extramural lesions within or adjacent to the tracheobronchial tree or gastrointestinal tract through the accessory channel of an ultrasound endoscope for Fine Needle Aspiration.

**Pentax** - Used to sample targeted submucosal and extramural lesions within or adjacent to the tracheobronchial tree through the accessory channel of an ultrasound endoscope for Fine Needle Aspiration.

This device is supplied sterile and is disposable - intended for single use only.

- High Definition dimpling pattern for exceptional needle visibility
- Naturally contoured handle offers needle control & stability
- Adjustable sheath length
- 5 cm needle extension presents a large diagnostic reach
- On scopes with metal luer hubs, the removable sheath adjuster allows direct connection to the scope
- On scopes with metal non-luer hubs, the sturdy scope adaptor design allows:
  - secure connection and procedural stability
  - suction and injection capability due to the airtight valve

Order Number	Reference Part Number	Needle Gage	Needle Extension cm	Sheath Fr	Dimple Pattern	Minimum Accessory Channel mm
G52011	*ECHO-HD-22-EBUS-O	22	5	4.1	High Definition	2.0
G53408	**ECHO-HD-22-EBUS-P	22	5	4.1	High Definition	2.0
G24062	*ECHO-HD-25-EBUS-O	25	5	4.1	High Definition	2.0
G24061	**ECHO-HD-25-EBUS-P	25	5	4.1	High Definition	2.0

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

\*Needle attaches to scope adapter on Olympus Endobronchial Ultrasound scope with non-luer metal hub.

\*\*Needle attaches to luer metal hub on Pentax Endobronchial Ultrasound scope, check for availability.

### Customer Service

EMEA: EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)

Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)

Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)

Belgium: +32 27001633, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)

Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)

Finland: +358 972519996, [fi.orders@cookmedical.com](mailto:fi.orders@cookmedical.com)

France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)

Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)

Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)

Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)

Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)

Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)

Norway: +47 23162968, [no.orders@cookmedical.com](mailto:no.orders@cookmedical.com)

Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)

Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)

Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)

Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)

Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)

United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

Americas: EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: [orders@cookmedical.com](mailto:orders@cookmedical.com)

Australia:

Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283

E-mail: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)



[www.cookmedical.com](http://www.cookmedical.com)