




Flexor[®] Ansel

GUIDING SHEATH



Track and treat with confidence.

REAL MANEUVERABILITY

- Developed for reliable kink resistance
- Hydrophilic coating smooths access
- Soft-tipped dilators enhance trackability

FLEXOR[®] ANSEL

- Many include .018" and .038" compatible dilators
- Straight, multipurpose, renal double and cobra curves available
- Features multiple lengths and sizes 4, 5, 6, 7, 8, 10 and 12 Fr

HIGH-FLEX ANSEL

- Stiffer shaft designed for improved pushability
- Softer, flexible dilator for unique trackability
- Straight and multipurpose curve options
- Available in 4, 5, 6, 7, 8, 10 and 12 Fr sizes in various lengths



Flexor[®] Ansel


GUIDING SHEATH

A pioneer in innovative access.

Where access is difficult, Flexor adapts. Where delivery is demanding, Flexor rises to the challenge.

Kink resistance

Patented coil-reinforced construction provides optimal flexibility and maximum resistance to kinking or compression.




Low friction

Large, low-friction PTFE-coated lumen facilitates device delivery.

Greater visibility

Radiopaque band, seamlessly incorporated into the sheath material, identifies precise location of the distal tip for positioning accuracy.



Enhanced trackability

Flexor Ansel has AQ[®] hydrophilic coating to dramatically reduce friction and enhance trackability.

Configurations

ANL0
Straight



ANL1
Multipurpose



ANL2
Renal Double



ANL3
Cobra




Flexor Ansel Guiding Sheath with Tuohy-Borst Side-Arm Adapter

Includes .018" and .038" compatible dilators



Tuohy-Borst Valve

Order Number	Reference Part Number	Introducer Fr	Sheath ID inch/mm	Accepts Wire Guide Diameter inch	Sheath Length cm	Tip Configuration
G29995	KSAW-5.0-18/38-55-RB-ANL1-HC	5.0	.074/1.9	.018/.038	55	multipurpose
G30120	KSAW-5.0-18/38-55-RB-ANL2-HC	5.0	.074/1.9	.018/.038	55	renal double
G30121	KSAW-5.0-18/38-55-RB-ANL3-HC	5.0	.074/1.9	.018/.038	55	cobra
G29996	KSAW-6.0-18/38-55-RB-ANL1-HC	6.0	.087/2.2	.018/.038	55	multipurpose
G29997	KSAW-6.0-18/38-55-RB-ANL2-HC	6.0	.087/2.2	.018/.038	55	renal double
G29999	KSAW-7.0-18/38-55-RB-ANL1-HC	7.0	.100/2.5	.018/.038	55	multipurpose
G30000	KSAW-7.0-18/38-55-RB-ANL2-HC	7.0	.100/2.5	.018/.038	55	renal double



Flexor Ansel Guiding Sheath with Check-Flo® Hemostasis Valve and High-Flex Dilator

Includes .035" compatible dilator



Check-Flo Valve

Order Number	Reference Part Number	Introducer Fr	Sheath ID inch/mm	Accepts Wire Guide Diameter inch	Sheath Length cm	Tip Configuration
G56229	KCFW-4.0-35-55-RB-HFANL0-HC	4.0	.059/1.5	.035	55	straight
G56234	KCFW-4.0-35-55-RB-HFANL1-HC	4.0	.059/1.5	.035	55	multipurpose
G56219	KCFW-4.0-35-70-RB-HFANL0-HC	4.0	.059/1.5	.035	70	straight
G56224	KCFW-4.0-35-70-RB-HFANL1-HC	4.0	.059/1.5	.035	70	multipurpose
G56230	KCFW-5.0-35-55-RB-HFANL0-HC	5.0	.074/1.9	.035	55	straight
G56235	KCFW-5.0-35-55-RB-HFANL1-HC	5.0	.074/1.9	.035	55	multipurpose
G56220	KCFW-5.0-35-70-RB-HFANL0-HC	5.0	.074/1.9	.035	70	straight
G56225	KCFW-5.0-35-70-RB-HFANL1-HC	5.0	.074/1.9	.035	70	multipurpose
G56231	KCFW-6.0-35-55-RB-HFANL0-HC	6.0	.087/2.2	.035	55	straight
G56236	KCFW-6.0-35-55-RB-HFANL1-HC	6.0	.087/2.2	.035	55	multipurpose
G56221	KCFW-6.0-35-70-RB-HFANL0-HC	6.0	.087/2.2	.035	70	straight
G56226	KCFW-6.0-35-70-RB-HFANL1-HC	6.0	.087/2.2	.035	70	multipurpose
G56232	KCFW-7.0-35-55-RB-HFANL0-HC	7.0	.100/2.5	.035	55	straight
G49043	KCFW-7.0-35-55-RB-HFANL1-HC	7.0	.100/2.5	.035	55	multipurpose
G56222	KCFW-7.0-35-70-RB-HFANL0-HC	7.0	.100/2.5	.035	70	straight
G56227	KCFW-7.0-35-70-RB-HFANL1-HC	7.0	.100/2.5	.035	70	multipurpose
G56233	KCFW-8.0-35-55-RB-HFANL0-HC	8.0	.113/2.9	.035	55	straight
G49045	KCFW-8.0-35-55-RB-HFANL1-HC	8.0	.113/2.9	.035	55	multipurpose
G56223	KCFW-8.0-35-70-RB-HFANL0-HC	8.0	.113/2.9	.035	70	straight
G56228	KCFW-8.0-35-70-RB-HFANL1-HC	8.0	.113/2.9	.035	70	multipurpose
G56278	KCFW-10.0-35-45-RB-HFANL0-HC	10.0	.134/3.4	.035	45	straight
G47700	KCFW-10.0-35-45-RB-HFANL1-HC	10.0	.134/3.4	.035	45	multipurpose
G56279	KCFW-12.0-35-45-RB-HFANL0-HC	12.0	.162/4.1	.035	45	straight
G47701	KCFW-12.0-35-45-RB-HFANL1-HC	12.0	.162/4.1	.035	45	multipurpose

Not all part numbers shown on this product sheet may be approved for sale in all jurisdictions. Consult with your local Cook representative or customer service center for details.



Flexor Ansel Guiding Sheath with Check-Flo Hemostasis Valve

Includes .018" and .035" or .038" compatible dilators



Check-Flo Valve

Order Number	Reference Part Number	Introducer Fr	Sheath ID inch/mm	Accepts Wire Guide Diameter inch	Sheath Length cm	Tip Configuration
G48186	KCFW-4.0-18/35-45-RB-ANL0-HC	4.0	.059/1.5	.018/.035	45	Straight
G48188	KCFW-4.0-18/35-90-RB-ANL0-HC	4.0	.059/1.5	.018/.035	90	Straight
G48191	KCFW-4.0-18/35-110-RB-ANL0-HC	4.0	.059/1.5	.018/.035	110	Straight
G44153	KCFW-5.0-18/38-45-RB-ANL0-HC	5.0	.074/1.9	.018/.038	45	Straight
G29981	KCFW-5.0-18/38-45-RB-ANL1-HC	5.0	.074/1.9	.018/.038	45	Multipurpose
G30117	KCFW-5.0-18/38-45-RB-ANL2-HC	5.0	.074/1.9	.018/.038	45	Renal Double
G48189	KCFW-5.0-18/38-90-RB-ANL0-HC	5.0	.074/1.9	.018/.038	90	Straight
G29988	KCFW-5.0-18/38-90-RB-ANL1-HC	5.0	.074/1.9	.018/.038	90	Multipurpose
G48192	KCFW-5.0-18/38-110-RB-ANL0-HC	5.0	.074/1.9	.018/.038	110	Straight
G44154	KCFW-6.0-18/38-45-RB-ANL0-HC	6.0	.087/2.2	.018/.038	45	Straight
G29982	KCFW-6.0-18/38-45-RB-ANL1-HC	6.0	.087/2.2	.018/.038	45	Multipurpose
G29983	KCFW-6.0-18/38-45-RB-ANL2-HC	6.0	.087/2.2	.018/.038	45	Renal Double
G29984	KCFW-6.0-18/38-45-RB-ANL3-HC	6.0	.087/2.2	.018/.038	45	Cobra
G29989	KCFW-6.0-18/38-90-RB-ANL1-HC	6.0	.087/2.2	.018/.038	90	Multipurpose
G44155	KCFW-7.0-18/38-45-RB-ANL0-HC	7.0	.100/2.5	.018/.038	45	Straight
G29985	KCFW-7.0-18/38-45-RB-ANL1-HC	7.0	.100/2.5	.018/.038	45	Multipurpose
G29986	KCFW-7.0-18/38-45-RB-ANL2-HC	7.0	.100/2.5	.018/.038	45	Renal Double
G49099	KCFW-8.0-18/38-45-RB-ANL0-HC	8.0	.113/2.9	.018/.038	45	Straight
G49100	KCFW-8.0-18/38-45-RB-ANL1-HC	8.0	.113/2.9	.018/.038	45	Multipurpose
G49101	KCFW-9.0-18/38-45-RB-ANL0-HC	9.0	.126/3.2	.018/.038	45	Straight
G49102	KCFW-9.0-18/38-45-RB-ANL1-HC	9.0	.126/3.2	.018/.038	45	Multipurpose

Customer Service Centers

AORTIC INTERVENTION
CRITICAL CARE
ENDOSCOPY
INTERVENTIONAL RADIOLOGY
LEAD MANAGEMENT
PERIPHERAL INTERVENTION
SURGERY
UROLOGY
WOMEN'S HEALTH

EMEA: EDI - www.cookmedical.com/edi.do
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Poland: +48 223060159, pl.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335,
 E-mail: orders@cookmedical.com

Australia:
 Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283,
 E-mail: cau.custserv@cookmedical.com



www.cookmedical.com