



Maxum[®]

REUSABLE FORCEPS



Get the surgical edge

When your procedure calls for maximum tissue sampling, the Maxum Reusable Forceps delivers. Surgical stainless steel jaws facilitate optimal sampling while the smooth coiled shaft enhances insertion and maneuverability through the endoscope. Whether performing a biopsy or a retrieval, maximize your results with Maxum.

www.cookmedical.com



Maxum[®]

REUSABLE FORCEPS

Used to obtain endoscopic mucosal tissue biopsies and/or for foreign body retrieval.

This device is supplied non-sterile and is reusable if device integrity is intact.

- Delivers maximum tissue sampling
- Dependable and reusable = cost effectiveness
- Sturdy, color-coded handle

ORDER NUMBER	GPN	CUP DIAMETER	FORCEP TYPE	SHAFT LENGTH	COLOR CODE	MINIMUM ACCESSORY CHANNEL
CBF-2.5-230	G22242	2.5 MM	NON-SPIKED	230 CM	RED	2.8 MM
CBF-2.5-230-S	G22243	2.5 MM	SPIKED	230 CM	RED	2.8 MM
GBF-1.8-160	G22244	1.8 MM	NON-SPIKED	160 CM	GREEN	2.0 MM
GBF-1.8-160-S	G22245	1.8 MM	SPIKED	160 CM	GREEN	2.0 MM
GBF-2.5-160	G22246	2.5 MM	NON-SPIKED	160 CM	GREEN	2.8 MM
GBF-2.5-160-S	G22247	2.5 MM	SPIKED	160 CM	GREEN	2.8 MM
RTF-1.8-160	G22253	1.8 MM	RAT TOOTH	160 CM	GREEN	2.0 MM
RTF-2.5-160	G22254	2.5 MM	RAT TOOTH	160 CM	GREEN	2.8 MM
RTF-2.5-230	G22255	2.5 MM	RAT TOOTH	230 CM	RED	2.8 MM



COOK ENDOSCOPY
 4900 Bethania Station Road, Winston-Salem, NC 27105 U.S.A.
 Phone: 336 744-0157, Toll Free: (USA) 800 457-4500, Fax: 336 744-5231
 Maxum is a registered trademark of Cook Endoscopy.
 Cook is a registered trademark of Cook Incorporated.

© 2009 Cook Medical

www.cookmedical.com