




OASIS[®]

ONE ACTION STENT INTRODUCTION SYSTEM



Streamline your plastic stenting procedures.

With its single-step delivery system, Cook's Oasis One Action Stent Introduction System effectively eliminates procedure steps, reducing exchange time. Whether you choose a conveniently preloaded Cotton-Leung Stent or ST-2 Tannenbaum, the Oasis delivers fast, efficient stent placement designed to streamline your plastic stenting procedures.

www.cookmedical.com



OASIS®

ONE ACTION STENT INTRODUCTION SYSTEM

Order Number	Reference Part Number	French Size Fr	Length cm
G25656	OACL-7-5	7	5
G25657	OACL-7-7	7	7
G25658	OACL-7-9	7	9
G25659	OACL-7-12	7	12
G25660	OACL-7-15	7	15
G25374	OACL-8.5-5	8.5	5
G25375	OACL-8.5-7	8.5	7
G25376	OACL-8.5-9	8.5	9
G25377	OACL-8.5-12	8.5	12
G25378	OACL-8.5-15	8.5	15
G25379	OACL-10-5	10	5
G25380	OACL-10-7	10	7
G25381	OACL-10-9	10	9
G25382	OACL-10-12	10	12
G25383	OACL-10-15	10	15
G25684	OACL-11.5-5	11.5	5
G25685	OACL-11.5-7	11.5	7
G25686	OACL-11.5-9	11.5	9
G25687	OACL-11.5-12	11.5	12
G25688	OACL-11.5-15	11.5	15

Set Includes: Introducer with four radiopaque bands 5 cm apart, with preloaded Cotton-Leung stent, 170 cm pushing catheter, 205 cm guiding catheter. Recommend Tracer Metro wire guide, sold separately. Length measured between flaps.

Order Number	Reference Part Number	French Size Fr	Length cm
G26846	OATS-8.5-5	8.5	5
G26848	OATS-8.5-7	8.5	7
G26849	OATS-8.5-9	8.5	9
G26853	OATS-8.5-12	8.5	12
G26854	OATS-8.5-15	8.5	15
G25658	OATS-10-5	10	5
G25659	OATS-10-7	10	7
G25660	OATS-10-9	10	9
G25663	OATS-10-12	10	12
G25664	OATS-10-15	10	15
G25836	OATS-11.5-5	11.5	5
G25837	OATS-11.5-7	11.5	7
G25838	OATS-11.5-9	11.5	9
G25842	OATS-11.5-12	11.5	12
G25843	OATS-11.5-15	11.5	15

Set Includes: Introducer with four radiopaque bands 5 cm apart, with preloaded ST-2 Soehendra Tannenbaum stent, 170 cm pushing catheter, 205 cm guiding catheter. Recommend Tracer Metro wire guide, sold separately. Length measured between flaps.

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

Customer Service Centers

- AORTIC INTERVENTION
- CRITICAL CARE
- ENDOSCOPY
- INTERVENTIONAL RADIOLOGY
- LEAD MANAGEMENT
- PERIPHERAL INTERVENTION
- SURGERY
- UROLOGY
- WOMEN'S HEALTH

EMEA: EDI - www.cookmedical.com/edi.do
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Poland: +48 223060159, pl.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335,
 E-mail: orders@cookmedical.com
Australia: ABN: 51 141259118
 Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283,
 E-mail: wacustserv@cookaust.com.au



www.cookmedical.com