

Access and orientation.




Tri-Tome pc[®]

TRIPLE LUMEN SPHINCTEROTOME

Conveniently fitted with a DomeTip[®] and dedicated wire guide and injection lumens, the Tri-Tome pc delivers the exceptional access and injection capability you need to allow you to focus on better patient outcomes. And, on our Protector[™] model, we've added insulation on the proximal portion of Tri-Tome's cutting wire to reduce the incidence of contact between the cutting wire and the endoscope.



www.cookmedical.com



Tri-Tome pc®

TRIPLE LUMEN SPHINCTEROTOME

Used for endoscopic cannulation of the ductal system and for sphincterotomy. Supplied sterile and is disposable - intended for single use only.

- Precurved DomeTip is designed with anatomy in mind for smooth, easy access
- Patented 3-D forming wire preserves optimal tip shape during shipping and storage, assuring superb orientation during the procedure
- Protector™ insulation on proximal portion of cutting wire decreases incidence of cutting wire contact with scope, reducing the risk of damage to the cutting wire and scope (available on TRI-25M-P)

Order Number	Reference Part Number	Tip	Cutting Wire Length mm	Cutting Wire Type	Wire Guide* Diameter inch	Color Connection
Tri-Tome pc® and Tri-Tome pc® Protector™						
G22583	TRI-20	DomeTip	20	Braided	.035	Purple
G22584	TRI-20M	DomeTip	20	Monofilament	.035	Purple
G22554	TRI-25	DomeTip	25	Braided	.035	Purple
G22555	TRI-25M	DomeTip	25	Monofilament	.035	Purple
G25502	TRI-25M-MET	DomeTip	25	Monofilament	.035	Purple
G26826	TRI-25M-SLT	DomeTip	25	Monofilament	.035	Purple
G22556	TRI-30	DomeTip	30	Braided	.035	Purple
G22557	TRI-30M	DomeTip	30	Monofilament	.035	Purple
G49197	TRI-25M-P	DomeTip	25	Monofilament	.035	Purple

Minimum accessory channel 2.8 mm. 7 FR catheter. *Recommend Tracer Metro® Wire, sold separately. Active cord available separately.

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

Customer Service Centers

EMEA: EDI - www.cookmedical.com/edi.do
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: orders@cookmedical.com

Australia:
 Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283
 E-mail: cau.custserv@cookmedical.com



www.cookmedical.com