



All the benefits of the Zilver® in a 6 Fr introducer.



Images courtesy of Dr. Jennifer Chennat
University of Chicago.

Zilver 635®

BILIARY SELF-EXPANDING STENT

Patient outcome is the only true way of judging a stent's quality. That's why the 6 Fr Zilver balances important stent traits - nonforeshortening, radial force, flexibility, durability and extraordinary patency rates when comparing the 10 mm Zilver to the 10 mm WallStent.¹ Being confident in a stent's performance allows you to focus on providing better results for your patient.



www.cookmedical.com



Zilver 635[®]

BILIARY SELF-EXPANDING STENT

Used in the palliation of malignant neoplasms in the biliary tree.

Supplied sterile and is disposable - intended for single use only.

Warning: The safety and effectiveness of this device for use in the vascular system have not been established.

- 6 Fr Flexor[®] introducer delivers maximum flexibility without kinking or compressing
- "Shelfless" tip designed for smooth withdrawal of introducer through deployed stent
- Reference markers on sheath and inner catheter for confident transpapillary placement
- Self-expanding, nonforeshortening nitinol stent provides maximum flexibility without kinking or compressing

Order Number	Reference Part Number	Stent Diameter mm	Stent Length cm	Introducer Fr	Introducer Length cm	Wire Guide* Diameter inch
G50616	ZILBS-635-6-4	6	4	6.0	200	.035
G50619	ZILBS-635-6-6	6	6	6.0	200	.035
G50622	ZILBS-635-6-8	6	8	6.0	200	.035
G50617	ZILBS-635-8-4	8	4	6.0	200	.035
G50620	ZILBS-635-8-6	8	6	6.0	200	.035
G50623	ZILBS-635-8-8	8	8	6.0	200	.035
G50618	ZILBS-635-10-4	10	4	6.0	200	.035
G50621	ZILBS-635-10-6	10	6	6.0	200	.035
G50624	ZILBS-635-10-8	10	8	6.0	200	.035

*Recommend Tracer Metro[®] wire guide, sold separately. Minimum accessory channel 2.8 mm.

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

¹ Loew BJ, Howell DA, Sanders MK, et al. Comparative performance of uncoated, self-expanding metal biliary stents of different designs in 2 diameters: final results of an international multicenter, randomized, controlled trial. *Gastrointestinal Endoscopy* 2009; 70(3): 445-453.

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: orders@cookmedical.com

Australia:
 Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283
 E-mail: cau.custserv@cookmedical.com

