

Use the next generation of a proven IVC filter.



Celect[®] Platinum

VENA CAVA FILTER

Improved visibility

- The four radiopaque platinum markers can verify placement accuracy immediately after you detach the filter.

Centering the filter during placement

- The struts on the legs of the Celect Platinum filter were designed to help you better center the filter during placement.
- The flexibility and trackability of Flex Tip on the femoral set helps you center the filter in the IVC.

Fixation

- Specially designed anchors and platinum markers on the filter's feet allow secure caval fixation.



www.cookmedical.com



Celect® Platinum

VENA CAVA FILTER

The product is intended for percutaneous placement via a vein for filtration of inferior vena cava (IVC) blood to prevent recurrent pulmonary embolism.

- Preadloaded on a filter introducer, the filter has four radiopaque platinum markers for accurate placement.
- The filter has a 7 Fr coaxial introducer system and a hydrophilic-coated 10 Fr predilator.
- The NavAlign™ multipurpose dilator includes radiopaque sizing bands and flushing sideports.

Order Number	Reference Part Number	Catheter Fr	Filter Max. Diameter mm	Filter Length mm	Catheter Length cm	Sheath Length cm
Femoral with NavAlign Delivery						
G34502	IGTCFS-65-1-FEM-CELECT-PT	7	30	49	79	65
Jugular with NavAlign Delivery						
G34309	IGTCFS-65-1-JUG-CELECT-PT	7	30	49	71	65
UniSet with NavAlign Delivery						
G34505	IGTCFS-65-1-UNI-CELECT-PT	7	30	49	71/79	65

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: orders@cookmedical.com

Australia:
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 E-mail: cau.custserv@cookmedical.com

