

# From cytology to histology: Advance your EBUS diagnostic procedures.



## EchoTip ProCore<sup>®</sup>

HD ENDOBRONCHIAL ULTRASOUND BIOPSY NEEDLE

At Cook, we were among the first to recognize the clinical impact of EBUS procedures. That's why we focused our expertise as market leaders in ultrasound needles to create a full range of EchoTip ProCore EBUS needles to meet your endobronchial ultrasound needs. Now, for the first time, you can achieve cytology and histology in your EBUS procedures, even in the challenging pulmonary hilar region. Our goal is simple: Give you the best histological information to help you deliver the best patient care.




[www.cookmedical.com](http://www.cookmedical.com)




Olympus



Pentax



25 gage



22 gage

## EchoTip ProCore®

HD ENDOBRONCHIAL ULTRASOUND BIOPSY NEEDLE

**Olympus** - Used with an ultrasound endoscope for fine needle biopsy, (FNB), of submucosal and extramural lesions within or adjacent to the tracheobronchial tree or gastrointestinal tract.

**Pentax** - Used with an ultrasound endoscope for fine needle biopsy, (FNB), of submucosal and extramural lesions within or adjacent to the tracheobronchial tree.

Supplied sterile and is disposable - intended for single use only.

- Featuring the core trap designed for receiving tissue in the delicate pulmonary area
- High Definition dimpling pattern for exceptional needle visibility
- Naturally contoured handle offers needle control & stability
- Adjustable sheath length
- 5 cm needle extension presents a large diagnostic reach
- On scopes with metal luer hubs, the removable sheath adjuster allows direct connection to the scope
- On scopes with metal non-luer hubs, the sturdy scope adaptor design allows:
  - secure connection and procedural stability
  - suction and injection capability due to the airtight valve

Order Number	Reference Part Number	Needle Size Gage	Sheath Size Fr	Adjustable Needle Extension cm	Stylet Tip	Minimum Accessory Channel mm
G34281	ECHO-HD-22-EBUS-O-C*	22	4.1	0-5	recessed ball tip	2
G34279	ECHO-HD-22-EBUS-P-C**	22	4.1	0-5	recessed ball tip	2
G34282	ECHO-HD-25-EBUS-O-C*	25	4.1	0-5	beveled	2
G34280	ECHO-HD-25-EBUS-P-C**	25	4.1	0-5	beveled	2

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

\*Needle attaches to scope adapter on Olympus Endobronchial Ultrasound scope with non-luer metal hub.

\*\*Needle attaches to luer metal hub on Pentax Endobronchial Ultrasound scope.

### Customer Service

**EMEA:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
 Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)  
 Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)  
 Belgium: +32 27001633, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)  
 Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)  
 Finland: +358 972519996, [fi.orders@cookmedical.com](mailto:fi.orders@cookmedical.com)  
 France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)  
 Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)  
 Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)  
 Norway: +47 23162968, [no.orders@cookmedical.com](mailto:no.orders@cookmedical.com)  
 Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)  
 Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)  
 Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

[www.cookmedical.com](http://www.cookmedical.com)

**Americas:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335  
 E-mail: [orders@cookmedical.com](mailto:orders@cookmedical.com)

**Australia:**  
 Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283  
 E-mail: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)

