Clip with instinctive control and precision.



Instinct[™] ENDOSCOPIC HEMOCLIP

During critical hemostasis procedures, you need a clip with unprecedented precision that can be simply and confidently controlled. With a handle and clip that function instinctively, the Instinct Hemoclip gives you full control of a 360-degree rotatable jaw for precise clip orientation. At 16 mm, the jaw is the widest on the market, is totally adjustable and features anchoring teeth so that you can securely grasp as much or as little tissue as you need. For your next clipping procedure, follow your instincts to enhanced patient care.





Instinct

ENDOSCOPIC HEMOCLIP

This device is used for endoscopic clip placement within the gastrointestinal tract for the purpose of endoscopic marking, hemostasis for mucosal/submucosal defects less than 3 cm in the upper GI tract, bleeding ulcers, arteries less than 2 mm, and polyps less than 1.5 cm in diameter in the GI tract. This device is not intended for the repair of GI tract lumenal perforations.

Supplied sterile and is disposable - intended for single use only.

- 360° jaw rotation provides complete control for precise clip placement
- **Wide** adjustable 16 mm opening span facilitates grasping variable amounts of tissue as desired
- Instinctive handle and clip design with a simple one-step clip detachment
- Opens and closes up to five times for optimal clip positioning
- Anchoring clip tips for secure tissue grasping
- MR Conditional per ASTM F2503 MR

Order Number	Reference Part Number	Length cm	Clip Opening mm	Minimum Accessory Channel mm	Box Qty
G18343*	INSC-7-230-S	230	16	2.8	10
G18344*	INSC-7-230-20-S	230	16	2.8	20
G18345*	INSC-7-230-40-S	230	16	2.8	40

*Check for availability.

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do Distributors: +353 61239240, ssc.distributors@cookmedical.com Austria: +43 179567121, oe.orders@cookmedical.com Belgium: +32 27001633, be.orders@cookmedical.com Denmark: +45 38487607, da.orders@cookmedical.com Finland: +358 972519996, fi.orders@cookmedical.com France: +33 171230269, fr.orders@cookmedical.com Germany: +49 6950072804, de.orders@cookmedical.com Hungary: +36 17779199, hu.orders@cookmedical.com Ireland: +353 61239252, ie.orders@cookmedical.com Italy: +39 0269682853, it.orders@cookmedical.com Netherlands: +31 202013367, nl.orders@cookmedical.com Norway: +47 23162968, no.orders@cookmedical.com Spain: +34 912702691, es.orders@cookmedical.com Sweden: +46 858769468, se.orders@cookmedical.com Switzerland - French: +41 448009609, fr.orders@cookmedical.com Switzerland - Italian: +41 448009609, it.orders@cookmedical.com Switzerland - German: +41 448009609, de.orders@cookmedical.com United Kingdom: +44 2073654183, uk.orders@cookmedical.com



www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335 E-mail: orders@cookmedical.com

Australia:

Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283 E-mail: cau.custserv@cookmedical.com