



# Quicksilver<sup>®</sup>

## BIPOLAR PROBE

### Smooth coagulation procedures

When you select the Quicksilver Bipolar Probe for your next coagulation procedure, expect smooth going. The Quicksilver's central lumen permits easy irrigation and the Hydromer<sup>®</sup> grafted tip reduces sticking for a smooth procedure – all at a competitive price.

[www.cookmedical.com](http://www.cookmedical.com)

AORTIC INTERVENTION

CARDIOLOGY

CRITICAL CARE

ENDOSCOPY

PERIPHERAL INTERVENTION

SURGERY

UROLOGY

WOMEN'S HEALTH



# Quicksilver®

## DISPOSABLE BIOPSY FORCEPS

Used in conjunction with a single or dual plug bipolar electro-surgical generator to endoscopically provide hemostasis throughout the gastrointestinal tract.

- Hydromer® grafted tip and spiral electrodes reduce sticking for effective coagulation
- Central lumen allows easy irrigation
- Extra long 350 cm length facilitates use with enteroscopes
- Available for use with single or dual plug bipolar generators for total compatibility
- Completely disposable, eliminating cleaning and reprocessing costs
- Full-body catheter for enhanced pushability

ORDER NUMBER	GPN	CATHETER SIZE	CATHETER LENGTH	GENERATOR PLUG TYPE
BCP-7A	G24925	7 FR	350 CM	DUAL
BCP-7B	G24926	7 FR	350 CM	SINGLE
BCP-10A	G24927	10 FR	350 CM	DUAL
BCP-10B	G24928	10 FR	350 CM	SINGLE



**COOK ENDOSCOPY**  
 4900 Bethania Station Road, Winston-Salem, NC 27105 U.S.A.  
 Phone: 336 744-0157, Toll Free: (USA) 800 457-4500, Fax: 336 744-5231  
 Quicksilver is a registered trademark of Cook Endoscopy.  
 Hydromer is a registered trademark of Hydromer, Inc.  
 Cook is a registered trademark of Cook Incorporated.

[www.cookmedical.com](http://www.cookmedical.com) © 2008 Cook Medical